PURA™ 4
PURA™ 10
PURA™ 14
PURA™ 22
PURA™ 30

Wasserbäder

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The above-named product is in compliance with the following harmonized standards and technical specifications:

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Water Baths

Operating instructions

Important: keep original operating manual for future use.

Table of contents

1 Use according to intended purpose	32
1.1 Description	32
2 Explanation of safety information.....	33
2.1 Explanation of other information	33
3 Operator's responsibility - safety instructions	34
3.1 Safety instructions for handling the device.....	34
3.2 Professional operation.....	35
3.3 Safety instructions for operation	35
4 Control and function elements.....	36
5 Safety instructions	38
6 Preparations	40
6.1 Installation	40
6.2 Mounting the bath cover (accessories).....	41
6.3 Bath fluid.....	41
6.4 Fill/empty.....	42
7 Initial operation	44
7.1 Instrument settings.....	45
8 Alarms, warnings / Possible causes.....	47
9 Technical Data	48
10 Cleaning / Repair of the device	52
10.1 Cleaning	52
10.2 Repair service	52
10.3 Send in for repair.....	52
11 Storage.....	53
12 Warranty	54
13 Waste disposal.....	55
13.1 Packaging.....	55
13.2 Unit	55
14 EC conformity	56

Congratulations.

You have made an excellent choice.

JULABO thanks you for the trust.

This operating manual is designed to familiarize you with the operation of our units and their possible applications. Please read the operating manual carefully.

Please call us if you have any questions about the operation of the unit or about the operating manual.

Contact:

JULABO GmbH
Gerhard-Juchheim-Strasse 1
77960 Seelbach/Germany
Tel. +49 7823 51-0
Fax +49 7823 2491
info.de@julabo.com
www.julabo.com

The JULABO quality management system

The standards for the development, production and distribution of temperature control devices for laboratory and industry use satisfy the requirements of ISO 9001 and ISO 14001. Registration certificate No. 01 100044846.

Unpacking and inspection

If the packaging is damaged or if you discover any concealed transport damage when you have unpacked the devices and the accessories, please notify the supplier in the form of a statement of damage.

NOTICE

The operating manual

- should be kept for future use.
- must be available to operating personnel at all times.

1 Use according to intended purpose

JULABO water baths have been designed for temperature application to specific fluids in a bath tank.

JULABO water baths are not conceived for direct temperature application to food and luxury articles or pharmaceutical and medico-technical products. Direct temperature application means: Unprotected contact of the object with the bath medium (bath fluid).

The unit is not for use in explosive atmosphere.

1.1 Description

PID1

- These water baths are operated using a (membrane) keypad which is protected from splashing water.
- The microprocessor technology allows the setpoint value to be displayed via the LED temperature display and set and stored via the (membrane) keypad.
- The PID temperature control automatically adjusts the heat supply to the desired setpoint temperature in the water bath.
- The water baths comply with the provisions of the European Directives applicable to the product.
- All PURA models have handles for transport.
- Models PURA 10, PURA 14, PURA 22 and PURA 30 are equipped with a drain screw for emptying.

2 Explanation of safety information

The operating manual contains warnings to increase safety when using the unit. The general warning sign, consisting of an equilateral triangle surrounding an exclamation sign and reproduced in various signal colors, is preceded by the signal words.

"Warning of a dangerous situation".

The significance of the danger is classified with a signal word. Read the instructions carefully and follow them.

DANGER

This signal word designates a danger with a high level of risk which, if it not prevented, will result in death or serious injury.

WARNING

This signal word designates a danger with a medium level of risk which, if it not prevented, may result in death or serious injury.

CAUTION

This signal word designates a danger with a low level of risk which, if it not prevented, may result in minor or moderate injury.

NOTICE

Designates a possibly harmful situation. If it is not prevented, the system or something near it may be damaged.

2.1 Explanation of other information

TIP

Your attention is drawn to something special by this.

Designates user tips and other useful information.

Dangers at second glance

Designates states which only occur after the start of an action and could have been prevented if the warning had been heeded.

Informative note

Provides additional information.

3 Operator's responsibility - safety instructions

Products manufactured by JULABO GmbH will operate safely if they are installed, operated and maintained as set out in the general safety regulations. This section explains the potential dangers which may occur when operating the water baths and sets out the main safety measures to prevent these dangers where possible.

The operator is responsible for the qualifications of its operating personnel.

Ensure that the personnel who operate the water baths are trained in the relevant work.

The operators must receive regular training about the dangers involved with their work and about action to prevent such dangers.

Ensure that everybody involved with the operation, maintenance and installation work have read and understood the safety information and the operating manual.

Please call us if you have any questions about using the design or about the operating manual!

3.1 Safety instructions for handling the device

- Avoid impacts on the housing, vibrations, damage to the operating keypad (keys and display) and heavy soiling.
- Ensure that the product is checked at regular intervals suitable for its frequency of use to ensure that it is in perfect condition.
- Check that instruction, warning, banning and safety marks are in good condition at regular intervals, but at least every 2 years.
- Ensure that the mains supply has a low impedance to prevent influencing other devices powered by the same supply.
- The device is designed for operation in a controlled electromagnetic environment. This means that in an environment of this nature, transmission equipment such as mobile phones should not be used in the immediate vicinity.
Other devices with components which are sensitive to magnetic fields may be influenced by magnetic radiation. We recommend that you maintain a minimum distance of 1 m.
- The ambient temperature must not exceed 40°C and not fall below 5°C.
- The relative humidity should not exceed 50% (40°C).
- Protect device from dirt.
- Protect the device from direct sunlight.

3.2 Professional operation

The water bath may only be configured, installed, maintained and repaired by trained personnel. The operator must be instructed by a specialist.

3.3 Safety instructions for operation

Only use the devices in well ventilated areas. The unit is not suitable for use in potentially explosive environments.

Do not pour flammable substances into the water bath.

Do not use **corrosive bath fluids** .

Thermal hazards due to the wide operating temperature range: hot, touchable parts and surfaces, superheated steam. Special care and attention is required.

When using hazardous substances or substances which may be dangerous, the operator must apply the enclosed safety symbols (1 + 2) on the front side of the device where they are clearly visible:

1		Warning of a danger zone. Attention! Observe documentation. (Operating manual, safety data sheet)
2a or		It is essential that you read the user information before switch on the device. Area of validity: EU
2b		It is essential that you read the user information before switch on the device. Area of validity: USA, NAFTA

4 Control and function elements

1 Recessed grip on both sides

2 Display field and control panel


2.1  LED temperature display, process value display / setpoint display (if decimal point flashes)
Alarm displays
Display of menu functions (parameters)

2.2  Control display - heating

2.3  Control display - alarm

2.4 °F °C Temperature units

2.5  Timer

2.6	h.m	Hours.minutes (timer display)
2.7		Edit buttons: Switch-over from actual value display to setpoint display. Set the temperature setpoint higher/lower. Press the key briefly for single steps. Keep key pressed for fast counting.
2.8		OK key 1. Switch the device on/off. 2. Save selected value / parameter.
2.9		MENU key
3		Drain screw (not available with PURA 4)
4		Mains switch
5		Mains connection
		Mains lead

5 Safety instructions

It is important to follow these safety instructions to prevent personal injury and property damage. These instructions apply in addition to the safety instructions at your workstation.

- It is essential that you read the user information before starting the device.
- Wear personal protective equipment (gloves, safety shoes, goggles).
- Only transport the device carefully and when it is empty. The interior of the unit can also be damaged by impacts or if it is dropped.
- For additional transport routes, transport the device on transport trolleys and secure it against falling down. Wear safety shoes.
- Do not loiter under the unit during transportation and operation.
- The water bath may only be operated with non-flammable fluids.
- Only operate the unit in rooms that are well-ventilated, dry and free of frost.
- Place the device on a flat surface and a mat made of non-flammable material.
- Do not start the device if it is damaged or leaking.
- Compare the mains voltage and frequency with the specifications on the type plate.
- Only connect the unit to a fused mains connection via a FI circuit breaker ($I_a=30$ mA).
- The mains plug is a safe isolator from the power supply network and must be freely accessible at all times.
- Check the mains cable for damage before initial operation of the device.
- Do not start the device if it has a damaged power cable.
- Escaping water vapor can reduce the operational safety of other devices in the vicinity of the water bath through condensation. This must be observed during installation and operation of the water bath.
- Parts of the bath cover can assume high surface temperatures at elevated working temperatures in continuous operation. If the bath cover is opened in this condition, there is an additional danger from the escaping hot water vapor.
- Caution when touching!
- Never use the device without bath fluid!
- Do not reach into the bath fluid.
- Do not drain the tempering medium when it is hot.

- Check the temperature of the bath fluid before draining it, for example by switching on the device briefly.
- Switch off the device and pull the plug before moving the device or carrying out service or repair work on the device.
- Have all service and repair work carried out by authorized specialists only.
- Switch off the device and disconnect it from the power supply before cleaning it.
- Empty the device completely before transporting it.

6 Preparations

6.1 Installation

Place the water bath on a horizontal and level surface.

WARNING

Risk of explosion

The device is not suitable for use in potentially explosive environments.

Hazards caused by the input material (samples)

The intended use of the water baths includes tempering samples in test tubes, Erlenmeyer flasks, etc., i.e. dipping them in.

We do not know what substances these vessels contain.

Many substances are:

- inflammable, flammable or explosive
- harmful
- polluting

in other words: dangerous

The user bears sole responsibility for handling these substances.

- Close sample containers properly if necessary.

NOTICE

Electrochemical oxidation, risk of corrosion on racks or samples made of non-ferrous metals.


- Avoid the use of insert racks and samples made of non-ferrous metal.
- Only use original JULABO insert racks.

6.2 Mounting the bath cover (accessories)

The bath cover is mounted with two hinges. The left side is shown.

5. Loosen the screws on the side of the device housing. Do not unscrew screws.

6. Push the hinges of the bath cover with the oblong holes over the screws.

7. Slide the hinges on the screws down to the stop.

8. Tighten the screws.
The bath cover is mounted.

6.3 Bath fluid

NOTICE

Risk of corrosion in the bath due to non-approved bath fluid!

Unsuitable, non-approved bath fluids can damage the water bath.

- Only water is permitted as a bath fluid.
- JULABO recommends a water mixture of 70% softened/descaled water and 30% tap water.
- Consult JULABO before using a bath fluid other than the recommended one.

- **Do not fill with ferrous water.** Rust formation also possible with stainless steel; damage to e.g. JULABO test tube racks.
- **Do not fill with chlorinated (bleached) water.** Pitting corrosion possible.
- **Do not fill with distilled and deionized water.** Corrosion in the bath possible, even with stainless steel.
- **Do not fill with hard water.** High lime content leads to calcification in the bath.
- **Do not use corrosive, or flammable bath fluids.**
- JULABO cannot accept any liability for damage caused by the use of an unsuitable bath fluid.

Water bath protection products

The water bath protection product "Aqua-Stabil" is recommended to combat algae, bacteria and fungus formation.

Order no.	Designation
8 940 006	6x 100 ml bottles
8 940 012	12x 100 ml bottles

6.4 Fill/empty

WARNING

Thermal hazards when opening the bath lid

Burns, scalds, hot steam, hot parts and surfaces which may be touched.

- Do not drain the tempering medium when it is hot.
- Allow the temperature of the bath fluid to cool down to room temperature before emptying.

Tip

The working height depends on the size and number of vessels to be inserted.

Fill the water bath only partially, insert the feed material and correct the filling level if necessary.

Tip:

Use bath lids to keep heat losses as low as possible, especially at working temperatures above 60°C.

Filling:

- Maximum filling height: 20 mm below the edge of the bath
Minimum filling height: 10 mm above the adjustable platform. Note the minimum marking on the platform.

Emptying:

7. Switch off the water bath at the mains switch.
 8. Disconnect the mains cable from the water bath.
 9. PURA 4 can now be distributed. The same applies to all other baths:
 10. Place a suitable receptacle for the bath fluid under the drain opening. Unscrew the drain screw.
 11. Empty the bath.
 12. After complete draining, close the drain opening with the drain screw.
- ✓ The bath is empty.

7 Initial operation

⚠ DANGER

Danger of electric shock if the housing is damaged!

If the device housing is damaged, electrically conductive parts may be touched.

- Do not operate or put into operation a device with damaged housing.

⚠ WARNING

Electrocution!

Risk of injury from electric power.

- Compare the mains voltage and frequency with the details on the model plate.
- Only connect the unit to a fused mains connection via a FI circuit breaker ($I_a=30$ mA).
- The device may only be connected to mains power outlets with a protective earth (PE)!
- The mains plug is a safe insulator of the power supply grid and must be freely accessible at all times.
- Do not start the device if it has a damaged power cable
- Check the mains cable regularly for signs of damage
- We cannot accept any liability for incorrectly connecting to the power supply

⚠ CAUTION

Danger from hot steam!

Danger of scalding when opening the bath cover

- Do not open bath cover at temperatures above 90°C with steam formation
- Allow device to cool down to room temperature
- Wear protective gloves and goggles
- When opening the bath cover, do not bend over the bath

NOTICE

Condensate formation due to steam escaping when bath cover is fitted!

Steam can escape through the vent holes and cause condensation to form in the direct vicinity of the water bath.

- Ensure sufficient distance between the water bath and adjacent devices and equipment
- Do not close vent holes

english

7.1 Instrument settings

	<p>Switch on the device -> Display OFF</p> <p>Press OK key -> display of actual temperature 28.5</p> <p>The temperature setting can be set when the device has been started or stopped.</p> <p>ⓘ The set value is saved so that it will be retained even after a power outage.</p>
<p>Setpoint adjustment</p> <p>Example:</p>	<p>Adjusting setpoints</p> <ol style="list-style-type: none"> 1. Press the edit key ▼ or ▲ briefly to switch from displaying the actual value to the setpoint. The decimal point will flash. 2. Change the value: Press the ▲ key to set a higher value. Press the ▼ key to set a lower value. Press the key briefly for single steps. Press and hold the key to adjust the values quickly. 3. Save the set value as the desired value by pressing the OK key. 4. The new setpoint will flash two times. A signal sounds at 0.5°C before the setpoint is reached.
<p>Timer setting</p>	<p>Setting the timer</p> <ol style="list-style-type: none"> 7. Press MENU key, "T" flashes OFF 8. Press key ▼ or ▲ key 0n 9. Press OK key run 10. Press the OK key 11. Use the ▼ or ▲ key to set h.m 12. Confirm with OK key, timer starts ("T" lights up).

	<p>After the setpoint is reached, the timer starts and the display alternates between the actual temperature and the remaining time. A second alarm signals at the end of the specified time.</p>
<p>Unit setting</p>	<p>Stopping the timer</p>
	<p>4. Press  key, "T" flashes On</p> <p>5. Press  or  key OFF</p> <p>6. Confirm with  key, timer stops 25.5 (actual temperature)</p>
<p>Set temperature units</p>	<p>4. Press key  2 x °C</p> <p>5. Press the  or  key °C or °F</p> <p>6. Confirm with  key, temperature values are displayed in the selected unit</p>

8 Alarms, warnings / Possible causes

For the malfunctions listed below, the unit will be permanently switched off.

The alarm indicator lights up and a continuous signal tone will also sound. The reason for the alarm or warning will be shown on the temperature display LED in coded form.

Warnings are displayed alternately with the actual value on the display.

The signal tone can be muted by pressing the OK key.

E 05

The cable for the working temperature sensor has been broken or short-circuited.

E 14

Bath temperature exceeds limit value of the temperature safety function. Cause: Fill level in the bath is too low and/or bath tank bottom is calcified. Remove the lime layer on the bath tank bottom. Observe the minimum fill quantity.

E 57

The bath temperature rises too quickly. Check the water level in the bath.

To cancel the alarm state

6. Switch off the unit at the mains switch.
7. Rectify the cause of the alarm.
8. Rectify the cause of the alarm or wait for approx. 4 seconds, depending on error type.
9. Switch on the device again at the main switch.
10. If the error occurs again, a remote diagnostic must be prepared.

JULABO Technical Service

Phone: +49 7823 5166

Fax: +49 7823 5166

Email: service.de@julabo.com

9 Technical Data

Water bath		PURA 4	PURA 10
Working temperature range ¹⁾	°C	18 ... 99.9	18 ... 99.9
Temperature stability ²⁾		±0.15	±0.15
Temperature setting		Digital	Digital
Temperature display		LED	LED
Display resolution	°C	0.1	0.1
Temperature Control		PID1	PID1
Heating capacity (at 230 V)	kW	0.5	1.2
Heating capacity (at 100-115 V)	kW	0.38 ... 0.5	0.9 ... 1.2
Bath opening (WxL)	cm	12x27	22x27
Bath depth	cm	14	14
Fill weight	liters	0.8 ... 4.5	1.4 ... 9.5
Overall dimensions (WxDxH) without bath cover	cm	21x35x22	31x35x22
Bath tank, basic material		Stainless steel	Stainless steel
Weight	kg	4.9	7.0
Ambient temperature	°C	5 ... 40	5 ... 40
Mains connection	V / Hz	230 ±10 % / 50 / 60	230 ±10 % / 50 / 60
Power consumption (at 230 V)	A	3	6
Mains connection	V / Hz	100-115 ±10% / 50/60	100-115 ±10% / 50/60
Current consumption (100 V / 115 V)	A	4 / 5	4 / 5

¹⁾ with counter-cooling/bath cover (accessories)

²⁾ with bath cover (accessories)

Water bath		PURA 14	PURA 22
Working temperature range ¹⁾	°C	18 ... 99.9	18 ... 99.9
Temperature stability ²⁾		±0.15	±0.15
Temperature setting		Digital	Digital
Temperature display		LED	LED
Display resolution	°C	0.1	0.1
Temperature Control		PID1	PID1
Heating capacity (at 230 V)	kW	1.8	2.0
Heating capacity (at 100-115 V)	kW	0.98 ... 1.3	0.98 ... 1.3
Bath opening (WxL)	cm	33x27	55x27
Usable bath depth	cm	14	18
Fill weight	liters	2 ... 14	3.4 ... 25.5
Overall dimensions (WxDxH) without bath cover	cm	42x35x22	64x35x26
Bath tank, basic material		Stainless steel	Stainless steel
Weight	kg	8.5	11.5
Ambient temperature	°C	5 ... 40	5 ... 40
Mains connection	V / Hz	230 ±10%/ 50/60	230 ±10%/ 50/60
Power consumption (at 230 V)	A	8	9
Mains connection	V / Hz	100-115 ±10%/ 50/60	230 ±10%/ 50/60
Current consumption (100 V / 115 V)	A	4 / 5	4 / 5

english

¹⁾ with counter-cooling/bath cover (accessories)

²⁾ with bath cover (accessories)

Water bath		PURA 30	
Working temperature range ¹⁾	°C	18 ... 99.9	
Temperature stability ²⁾		±0.15	
Temperature setting		Digital	
Temperature display		LED	
Resolution	°C	0.1	
Temperature Control		PID1	
Heating capacity (at 230 V)	kW	2.0	
Heating capacity (at 100-115 V)	kW	0,98 ... 1,3	
Bath opening (WxL)	cm	77x27	
Usable bath depth	cm	18	
Fill weight	liters	4.8 ... 36	
Overall dimensions (WxDxH) without bath cover	cm	86x35x26	
Bath tank, basic material		Stainless steel	
Weight	kg	14.5	
Ambient temperature	°C	5 ... 40	
Mains connection 230 V / 50 / 60 Hz	V / Hz	230 ±10%/ 50/60	
Power consumption (230 V)	A	9	
Mains connection	V / Hz	100-115 ±10%/ 50/60	
Current consumption (100 V / 115 V)	A	10 / 12	

¹⁾ with counter-cooling/bath cover (accessories)

²⁾ with bath cover (accessories)

All data refers to:

- Rated voltage and rated frequency
- Operating temperature 70°C
- Ambient temperature 20°C
- Water bath fluid

Safety precautions to IEC 61010-2-010:
 Safety temperature / overheating protection - fixed 130°C
 Classification to DIN 12876-1 Class I
 Alarm Optical and acoustic (permanent)

Ambient conditions to IEC 61010-1:

For interior use only.
 Altitude up to 200 m - normal zero.
 Ambient temperature: +5 ... +40 °C
 Humidity
 maximum relative humidity 80% for temperatures up to 31°C,
 decreasing linearly up to 50% relative humidity at 40°C
 Voltage fluctuations of ±10% are permitted.
 Protection type to EN 60 529: IP 21
 The device complies with safety class I
 Overvoltage category II
 Contamination level 2

EMC requirements according to EN61326-1:

The device is an ISM device of group 1 per CISPR 11 (uses HF for internal purposes) and is classified in class A (industrial and commercial sector).

NOTICE

- Devices of class A are intended for the use in an industrial electromagnetic environment.
- When operating in other electromagnetic environments, the electromagnetic compatibility with other devices may be impacted.

10 Cleaning / Repair of the device

CAUTION

**Danger of injury during maintenance, repair and transport
Danger from mains voltage.**

- Switch off the device and disconnect it from the power supply before undertaking the cleaning work.
- It is imperative that no moisture penetrate the inside of the circulator.
- Do not use alcohol- or solvent-based cleaning agents.
- Have all service and repair work carried out by authorized electricians only.

10.1 Cleaning

Regularly remove limescale deposits on the bath tank bottom with a limescale remover. Clean the functional parts of the device with a soft cloth and mild soapy water.

JULABO water baths are designed for continuous operation in normal conditions. No regular maintenance work is required.

The bath tank should only be filled with suitable bath fluid. In the event of contamination, the bath fluid must be replaced from time to time.

If Aqua Stabil is used, blue residues on the inside of the bath can be removed with a cleaning agent based on anionic surfactants (e.g., dish detergent). For this purpose, fill the bath with a diluted cleaning agent solution and operate it for several hour at 60 °C.

10.2 Repair service

Before calling out a service technician or sending back a JULABO device for repair, we recommend that you contact our Technical Service Department.

JULABO Technical Service

Phone: +49 7823 5166

Fax: +49 7823 5199

E-mail: service.de@julabo.com

10.3 Send in for repair

In the event that you send a device to JULABO:

- Clean the device and, if necessary, detoxify it appropriately to avoid the risk of harming the service personnel.

- Ensure that you pack the device carefully and correctly.
- It is essential that you enclose a brief description of the error.
- Before returning the device, please complete an online return form at <http://www.julabo.com/com/support/rma>.
- JULABO cannot accept any liability for damage caused by incorrect packaging.
- JULABO reserves the right to complete any necessary modifications of a technical type which help to ensure perfect function during the repair process in order to improve the product.

11 Storage

Devices that are not to be reused must be stored in a dry, place, protected from dust and frost, after cleaning. The system components must be fully emptied and carefully dried, for example using compressed air. Do not store in an aggressive atmosphere.

12 Warranty

JULABO provides a warranty that the device will function perfectly as long as it is connected and handled correctly and as described in the operating manual.

**The warranty period is
one year.**

Free extended warranty

2 Years Warranty

1Plus Warranty

Registration free of charge on www.julabo.com

The 1PLUS guarantee gives the user a free extended warranty to 24 months, limit to a maximum of 10,000 hours of service.

The requirement for this is that the user registers the device on the JULABO website at www.julabo.com, quoting its serial number. The JULABO GmbH invoice date is used to define the warranty period.

In the event of a justified complaint the warranty is limited to refinishing or free repair or replacement at our discretion. Defective parts will be repaired or replaced free of charge if, in the event of a fault or a defect, there is evidence that it was due to a material or workmanship flaw.

All other compensation claims will be rejected.

13 Waste disposal

13.1 Packaging

Packaging materials must be disposed of as prescribed by the current local regulations.

13.2 Unit

In the European Economic Area (EEA) the disposal of waste equipment is regulated in the "**Directive of the European Parliament and of the Council on Waste Electrical and Electronic Equipment (WEEE)**". The current official journal on this matter is available on the European Parliament's homepage.

The symbol for the separate collection of electrical and electronic equipment is a crossed-out trash can.

Disposal with household waste (unsorted waste) or similar collections of municipal waste is not permitted!

Contact an authorized waste disposal contractor in your country.

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JULABO USA, Inc.