# DEXIS™

User Manual  
Manuel d'utilisatin  
Manuale dell'utente  
Benutzerhandbuch  
Manual del usuario  
Gebruikershandleiding  
Manual do Utilizador

# DEXIS™

User Manual

---

# Table of Contents

## Chapter 1 - Introduction

Essential Performance .....	1-1
Indications for Use .....	1-2
Product Description .....	1-2
Proficiency with the DEXIS Sensor .....	1-3
Conventions Used in the Manual .....	1-4
Unpacking the System Components .....	1-4
Contents of Sensor Box .....	1-5
System Components .....	1-6
Digital Intraoral Sensor .....	1-6
Device Labeling .....	1-7
Cables .....	1-8
Sensor Garage .....	1-8
Documentation .....	1-8
USB Software Kit .....	1-8

## Chapter 2 - Safety and Disposal Procedures

Sensor Inspection .....	2-2
Electrical Safety .....	2-2
Mains Isolation .....	2-2
X-ray Protection .....	2-3
Prevention of Cross-Contamination .....	2-3
Product Disposal .....	2-3
Prevention of Environmental Contamination .....	2-3
Prevention Against Cybersecurity Threats .....	2-4

## Chapter 3 - System Configuration

Personal Computer Requirements .....	3-1
--------------------------------------	-----

## Chapter 4 - Installation

Installation Steps .....	4-1
Installation for Users without USB Software Kit .....	4-3
Gesture and Dashboard Feature Support .....	4-3
Dashboard Icons .....	4-5

Frequently Asked Questions for Installation .....4-6  
Mounting the Sensor Garage .....4-7

## Chapter 5 - Use

Acquisition of Radiographic Images .....5-1  
Portability .....5-2  
Hygiene .....5-2  
    Recommended Disinfecting Solutions .....5-3  
Maintenance .....5-3  
Image Quality Assurance .....5-3  
    X-ray Dose / Exposure Time Settings .....5-4  
    X-ray Sharpness / Contrast .....5-4  
    Display Image .....5-4

## Chapter 6 - Specifications and Standards

Sensor Specifications .....6-1  
Compatibility with Radiographic Generators .....6-2

## Chapter 7 - Troubleshooting

System Does Not Acquire X-ray Images .....7-1

## Appendix A - Accessories

## Appendix B - EMC Information

## Appendix C - Product Symbols

# Chapter 1 Introduction

This User Manual contains original instructions by Dental Imaging Technologies Corporation for safe set-up, use and maintenance of the DEXIS sensor. The user must be able to read and understand the written language of the manual.

It also contains technical specifications of the system and basic information on how the system works.

Please read this manual carefully before starting to use the device, paying particular attention to warnings, especially safety warnings. Operators are strongly urged to comply with the current recommendations of the International Commission on Radiological Protection and, in the United States, the US National Council for Radiological Protection.

## Essential Performance

### **WARNING**

---

In the event of a sensor malfunction or a noticeable change in image quality, it is recommended to immediately cease the use of the device on the patient and call Technical Support at 1-888-883-3947.

---

### **WARNING**

---

No modification of this equipment is allowed.

---

The essential performance of the DEXIS sensor system comprises:

- The ability of the sensor to capture X-ray images suitable for recognition of normal anatomical structures, dental pathologies, and abnormal conditions, where inadequate images may result in mis-diagnosis, subjecting the patient to incorrect or unnecessary dental procedures that would present an unacceptable risk to the patient.
- The ability of the sensor holder to align the sensor with the X-ray source for imaging the desired anatomy.
- The ability of the single use sensor sheaths to isolate the sensor from the patient, where cross-contamination could present an unacceptable risk to the patient.

## Indications for Use

The DEXIS sensor is a USB-driven digital sensor which is intended to acquire dental intraoral radiographic images. The DEXIS sensor shall be operated by healthcare professionals, who are educated and competent to perform the acquisition of dental intra-oral radiographs. The DEXIS sensor can be used either in combination with special positioning devices to facilitate positioning and alignment with the X-ray beam or it may also be positioned by hand with the assistance of the patient.

---

**WARNING**

Some components of the DEXIS Sensor may have the presence of Nickel greater than 0.1% by article weight. Contact may result in an allergic reaction due to Nickel sensitivity. Use of a sheath that covers the sensor positioning system may reduce the risk of exposure to the patient. Use of gloves while handling the alignment bar may reduce the risk of exposure to the operator.

Exposure to nickel may possibly cause cancer. Use of a sheath that covers the sensor positioning system may reduce the risk of exposure to the patient. Use of gloves while handling the alignment bar may reduce the risk of exposure to the operator.

---

---

**WARNING**

Take the necessary steps to protect yourself from radiation. For proper operator positioning, refer to the 'Instructions for Use' of your intraoral X-ray equipment.

---

---

**WARNING**

Under no circumstances should the dental professional hold the sensor by hand during the X-ray exposure.

---

## Product Description

The DEXIS Sensor is an indirect converting x-ray detector, e.g. incident X-rays are converted by a scintillating material into (visible) light, this light is coupled optically to a light detection imager based on CMOS technology.

The design of the sensor assembly supports the automatic detection of the incident x-rays to generate digital images for dental intra oral applications.

The DEXIS Sensor supports USB2.0 and USB3.x connectivity to personal computers using a dedicated electronic assembly and a sensor software driver.

## Proficiency with the DEXIS Sensor

Please become acquainted with your DEXIS sensor. The sensor is your tool. Learn to use it well and it will become an effective aid in your diagnosis, treatment and treatment planning process. Visual impairment is permissible so long as the user can resolve the necessary details of the accompanying documentation, user interface and image data. Use of corrective measures such as glasses or contact lenses is permissible. Audio impairment is permissible so long as the user can discern audio indications emitted by other associated equipment, such as the exposure buzzer that sounds while X-rays are being emitted, or, in the case of deafness, so long as such audio indications are supplemented by corresponding visual indications. Use of corrective measures such as hearing aids is permitted.

As with all new clinical tools, it is necessary to invest a certain amount of time for practice in order to become proficient with the sensor.

For new DEXIS users, we highly recommend that you register for in-office training and that you and all of your clinical team set aside time to learn the system together. This manual is not intended to serve as a substitute for this training. However, as preparation for the training, we do encourage you to read the manual thoroughly and practice with your equipment safely and carefully for at least one week prior to your training date.

Please note that this manual assumes that new users possess basic computer skills and an understanding of the Windows<sup>®</sup> operating system. Absent this experience, we highly recommend that you obtain these skills through a computer course, video, or textbook. Your imaging representative may be able to suggest (although not endorse) one or more of these computer learning resources.

---

## Conventions Used in the Manual

The following conventions are used to bring the operator's attention to important information:

	<b>Warning</b> Alerts the operator that failure to follow the procedure could cause bodily injury or death.
	<b>Caution</b> Alerts the operator that failure to follow the procedure could cause damage to the equipment or cause loss of data.
<b>Important:</b>	Provides advice for the operator regarding use of the device or a process.
<b>NOTE:</b>	Highlights important or unusual points.

## Unpacking the System Components


The DEXIS System is carefully inspected and packaged prior to shipment. If the DEXIS system was shipped to you, please remove the contents of the shipping container and be sure to identify and directly locate each of the system components shown below.

**NOTE:** Report any damaged components to the shipping company and any missing components to your dealer within 24 hours of receiving the shipment.


Prior to use, clean and disinfect sensor, cable and holders with CaviWipes.

---

## Contents of Sensor Box

Top Layer


Bottom Layer

## System Components

### Digital Intraoral Sensor





The DEXIS Sensors are USB-driven digital intraoral sensors based on a CMOS technology that has been specifically designed for dental applications. The DEXIS Sensors are available in three formats:

- DEXIS Titanium Sensor (active area 26x33mm)
- DEXIS IXS Size 1 (active area 20x30mm)
- DEXIS IXS Size 2 (active area 26x36mm)

# Device Labeling

## Sensor ID Labels (typical)

## Cables

One meter (3 feet) USB extender cable.

## Sensor Garage

Designed to provide a safe and secure storage solution for DEXIS Sensors.

## Documentation

1. DEXIS Sensor User Manual
2. DEXIS IXS Software Update Driver

## USB Software Kit

Contains update software drivers.

# Chapter 2 Safety and Disposal Procedures

The device must be installed and used in accordance with the safety regulations and instructions for use supplied in this User Manual, for the purposes and applications for which it is intended.

**Important:** Any serious incident that has occurred in relation to the use of the DEXIS Sensor should be reported to Technical Support at 1-888-883-3947.

All external surfaces of the sensor, sensor cable, sensor holders and the sensor sheaths are considered to be applied parts and are safe for normal or accidental patient contact during use.

The DEXIS Sensor has no serviceable parts and does not require calibration. Do not open the device to service it. All aspects of the DEXIS Sensor that are meant to be attended to by the operator are accessible without opening the internal components of the device. If there is a service problem, contact a qualified dealer service representative or Dental Imaging Technologies Corporation Support.

**Important:** Portable/mobile radio frequency communications equipment can affect the function of the DEXIS Sensor as well as any other electronic medical equipment.

DEXIS Sensor is a USB compliant device and shall be used with USB compliant cables suitable for high speed/USB 2.0 cables. Such cables are either marked "USB 2.0" or "USB Hi-Speed."

USB certified hubs or active booster cables can be used to extend the distance to the USB host/computer. The length of the cable connection to the hub or between hubs shall not exceed 5m.

**NOTE:** 3 ft. Booster Cable (1.014.1145) is available for purchase and sold separately.

## **CAUTION**

---

Using non-USB compliant cables or hubs, or exceeding the maximum count of USB hub devices for extending the distance, can degrade the immunity of the DEXIS Sensor to electromagnetic fields or increase the emission of electromagnetic fields from the sensor.

---

## Sensor Inspection

Always inspect the sensor, cable, and positioning devices for physical damage prior to every use.

### **WARNING**

---

Remove sensor from service if damage to the cable or housing is observed, otherwise exposure to elevated surface temperature or electrical energy may result.

---

## Electrical Safety

The product must be used only in rooms or areas which comply with all laws and regulations applicable to electrical safety in medical premises, where provisions for safety ground connections are present.

All X-ray equipment for Dental Intra-oral Radiography used with the DEXIS Sensor must conform to IEC 60601-2-65.

The DEXIS Sensor conforms to safety standard IEC 60601-1.

The DEXIS Sensor is not suitable to be operated in oxygen rich and/or explosive environments.

All IT components electrically connected to the DEXIS Sensor must conform to IEC 60950-1.

Normally, the IT components are placed OUTSIDE the patient environment. IT components placed INSIDE the patient environment, due to customer site requirements, must also conform to IEC 60601-1.

IEC 60601-1 defines the “Patient environment” as “any volume in which intentional or unintentional contact can occur between a Patient and parts of the ME Equipment or ME System or between a Patient and other persons touching parts of the ME Equipment or ME System.”

## Mains Isolation

Disconnection from the supply mains occurs at the input to the computer. The sensor can also be disconnected from the computer.

## X-ray Protection

The rules of dental radiography still apply to digital X-ray systems. Please continue to use protection for your patients. As a clinician, clear the immediate area when exposing the sensor.

### **CAUTION**

---

US Federal law restricts this device to sale by or on the order of a dentist or other licensed practitioner.

---

## Prevention of Cross-Contamination

To help prevent cross-contamination between patients, place a new hygienic barrier on the sensor for each new patient. The hygienic barrier must cover the sensor and at least 3-4 in. (7-10 cm) of the cable.

## Product Disposal

Dispose of sheaths and other consumables following the normal dental office procedure for biomedical waste. Improper disposal of biomedical waste can lead to the spread of illness and disease.

■ Please contact your dealer or supplier for further information about product disposal at the end of the product's lifetime.

## Prevention of Environmental Contamination

Dispose of sheaths and other consumables following the normal dental office procedure for biomedical waste. Improper disposal of biomedical waste can lead to the spread of illness of disease.

## Prevention Against Cybersecurity Threats

Protecting your practice against cybersecurity threats is shared responsibility between us as the manufacturer and you as the health care provider. DEXIS has taken precautions to ensure that your sensor, as shipped from the factory, is protected against such threats. Additionally, during on-site repairs, your DEXIS Sensor is only accessible by trained DEXIS Customer Support Representatives. However, we **cannot** protect your network system.

We recommend having a **qualified** IT security professional **hired by you** to check your network system for proper protection against viruses, malware, and intrusions (e.g. anti-virus software, and/or use of firewall).

### **WARNING**

---

Failure to maintain cyber security can result in compromised device functionality, loss of data (medical or personal) availability or integrity, or exposure of other connected devices or networks to security threats.

---

# Chapter 3 System Configuration

## Personal Computer Requirements

	Recommended
Operating System	Microsoft Windows 10 Professional 64
Processor	Core I5
Memory	8 GB
Hard Disk	512 GB
Display Setting	2560 x 1440 (QHD), 1920 x 1080 (minimum)
Video Memory	128+ MB
Ports	USB 2.0 / USB 3.x
Standards	In compliance with the UL/IEC/EN 60950-1 standard

**NOTE:** The DEXIS Sensor is a USB certified device and shall be used with USB compliant cables suitable for high speed/USB 2.0 components. Certified USB cable extenders are available from DEXIS. Powered USB active hubs can also be used to extend the distance to the USB host/computer. The total length of all cabling, including the DEXIS Sensor shall not exceed 5m, unless a powered USB active hub is used.

### **CAUTION**

Using non-USB compliant cables or hubs, or exceeding the maximum count of USB hub devices for extending the distance, can degrade the immunity of the DEXIS Sensor to electromagnetic fields or increase the emission of electromagnetic fields from the sensor.

### **Hard Disk**

The choice of hard disk capacity depends on the number and size of images to be stored. Intraoral images vary in size from 4.7 MB (uncompressed TIFF images with the Size 2 sensor) to 100 KB (JPEG/JFIF with 50% compression). To guarantee data security, however, it is preferable to allow for duplication to provide back-ups of files containing diagnostic images and information.

Example calculation: 10 patients per day x 20 images per patient x 250 workdays per year x 5 MB image= 250 GB per year.

## Backup

It is highly recommended to maintain a spare backup copy to permit data recovery. The database of patients and images should be copied frequently (e.g. once a week) onto a removable mass storage device (removable hard drive, USB, etc.). You may use the back-up functions offered by Microsoft® software or directly copy files of data and images onto a removable device.

## Supported Imaging Software

- VixWin Platinum 3.6.3 or higher
- DEXIS 9.5.0 or higher
- DEXIS 10.1.6 or higher
- DTX Studio Clinic 1.6 or higher
- CLINIVIEW 11.5.2 and DEXIS Driver 20.1 or higher
- GxPicture 4.0.1 or higher
- Instrumentarium Dental TWAIN 7.5 and DEXIS Driver 20.1 or higher

## Application Software

Some diagnostic imaging software, practice management software or other third-party software may have different minimal system requirements. Refer to the requirements of the specific software in its User Manual.

**Important:** If the DEXIS System is used with imaging software other than those listed above, then the Intended Use of that software needs to include dental intra oral X-ray.

# Chapter 4 Installation

## Installation Steps

**IMPORTANT:** Must be completed before using your sensor.

---

### CAUTION

Some dental imaging software uses a 3rd party SDK to integrate with DEXIS. Contact your software manufacturer or DEXIS Technical Support for software and installation instructions.

---

Your DEXIS Sensor comes with the Software Manager on the USB Software Kit. Simply insert this thumb drive and click the **Software Manager** shortcut to get started.

You will need a valid Internet connection. If your workstation is not connected to the Internet, you will need to follow the steps below to download the software and copy it to your workstation.

If you no longer have access to your USB Software Kit, you will need to follow the steps below to download the software and copy it to your workstation.

1. Log into your computer with an Administrator account.
2. Open your supported Imaging Software and verify that you have the minimum version installed. If installing a new Imaging Software for the first time, it is recommended to complete the installation prior to installing your sensor drivers.
3. Your DEXIS Sensor comes with the Software Manager on the USB Software Kit. Simply insert this thumb drive and click the Software Manager shortcut to get started.

---

### CAUTION

You will need a valid Internet connection to perform the installation. If your workstation is not connected to the Internet, you will need to contact DEXIS Technical Support.


If you no longer have access to your USB Software Kit, you will need to follow the steps below to download the software and copy it to your workstation.

---


- From the USB drive, right click on DEXIS\_IO\_Sensor\_Driver\_Setup\_Launcher.exe and select Run As Administrator.

The Software Manager is ready to install the latest version of your Imaging Software, the DEXIS Legacy IO Sensor Drivers (if not previously installed) and DEXIS IO Sensor Drivers.


- The Software Manager begins installing the software.

- Accept End User License Agreement.

- Select the Typical option to install the standard feature set.
- Verify that your Imaging Software and DEXIS IO Sensor Drivers were successfully installed by viewing the Software Manager.
- Reboot your computer to complete the software installation for your new DEXIS sensor

**NOTE:** The calibration of the sensor is initiated when a patient is launched in DEXIS 9 and when IO acquisition is launched in DEXIS 10. This is not needed for VixWin/GxPicture, as calibration begins automatically when the sensor is connected.

---

## Installation for Users without USB Software Kit

1. Log into your computer with an Administrator account.

Open your supported Imaging Software and verify that you have the minimum version installed.

2. Navigate to <https://dexis.com/en-us/download-center> and download the latest version of your Imaging acquisition software, according to the version that you are currently using.
3. Download the latest version of the DEXIS IO Sensor Drivers by searching for DEXIS Titanium. The file name for the driver is DEXIS\_IO\_Sensor\_Drivers.exe.

Open your Imaging acquisition software to install your update.


**NOTE:** You will NOT need to uninstall your previous version of your Imaging acquisition software in advance.

4. Right click the DEXIS\_IO\_Sensor\_Drivers.exe and select Run as Administrator.
5. Select the Typical option to install the standard feature set.
6. Reboot your computer to complete the software installation for your new DEXIS sensor.


**NOTE:** The calibration of the sensor is initiated when connected to your imaging software.

## Gesture and Dashboard Feature Support


1. Select Custom for Dashboard and/or Gesture support.


2. Select the desired program drivers.

3. Select DEXIS Titanium Gesture services (Titanium only).


4. The selected programs are installed.


5. After the InstallShield Wizard completes the installation of the sensor drivers, click Finish.

**NOTE:** If your Imaging Software's status is listed as Modified, the software was successfully installed but you can still install additional purchased software add-ons.

- Click Yes to restart your system to complete the software installation for your new DEXIS sensor.

## Dashboard Icons

The sensor dashboard provides status information. After installation and reboot are complete, a sensor with a white background and a red "X" is displayed on the right side of the Windows Task Bar. The icon is indicating that the sensor is not plugged into the USB port. This icon changes based on the status of your device.

**Red X** - indicates either:

- No DEXIS Sensors are connected to computer
- Or one of the connected sensors is in an error state.

**Green** - indicates that all of the connected sensors are ready to operate.


**Yellow** - indicates that at least one of the connected sensors is connecting or downloading calibration files.

**Blue** - indicates that one of the connected sensors is acquiring an image.


- Connect your sensor to your computer's USB port to start the automatic download of the calibration files.

**NOTE:** Right click this icon to open the *Intra-Oral Sensor Dashboard*.

- The DEXIS Sensor has the calibration files embedded on the device. Installation of the calibration files typically takes around three minutes to complete. You may open the *Intra-Oral Sensor Dashboard* to see the software load progress.

**NOTE:** DO NOT unplug the sensor while the calibration files are being downloaded.

- Once the download of calibration files are complete, the Sensor Dashboard icon background turns green.
- The Intra-Oral Sensor Dashboard also shows the serial number and state of the sensor. This screen also displays driver and software versions.

The sensor is now ready for use in Imaging Software.

## Frequently Asked Questions for Installation

**Question:** Do I need Windows Administrator rights and privileges to install the software and driver?

**Answer:** Yes. If a Windows user account does not have Administrator rights, the software and drivers may not install properly for the DEXIS Sensor

**Question:** Do I need to set exceptions/exclusions in my anti-virus/security programs?

**Answer:** Yes. Not doing so may cause the software and drivers to not operate correctly. The following folders need to be added for exceptions/ exclusions depending on which drivers and imaging software you are installing.

### DEXIS Driver

- C:\ProgramData\DEXIS
- C:\Program Files (x86)\DEXIS

### DEXIS 9

- C:\DEXIS

### VixWin

- C:\VixWin

---


**Question:** Why am I receiving an Error 72 message?

**Answer:** Inconsistent USB power or data from a USB port can interrupt the PolarisAgent Service. Using a powered USB 2.0/3.0 hub could resolve these issues.

## Mounting the Sensor Garage

Install the Sensor Garage in a place that is convenient and does not exceed the maximum length of the USB cable, which is 10 ft. (3 m). The sensor garage should also be positioned so that the cable from the sensor does not create a tripping hazard.

The Sensor Garage comes with two easy to install options. Use with the supplied double-sided tape backpaper or mount with drywall anchors and screws (not supplied). Remove the double-sided tape and firmly press the Sensor Garage against the mounting surface. If using screw to mount, use the two indicator notches on the top edge of the garage for proper spacing and alignment.

# Chapter 5 Use

## Acquisition of Radiographic Images

Turn on the PC that has the DEXIS Sensor installed and launch the imaging software, such as DTX or DEXIS Imaging Suite (please refer to the software manual for information on the software program).

1. Carefully plug the USB end of the sensor cable into the PC USB port. Do not use force, ensure that the plug is properly oriented with the PC port. Try rotating the USB plug if it does not easily insert into port.

### CAUTION

---

Use extreme care when mounting the sensor inside a cabinet. Pinching/crushing the sensor cord or USB Extension Cable can damage the cord and cause image derogation or failure.

---

2. Set the required technical settings (exposure time, etc.) on the X-ray generator (refer to X-ray Dose / Exposure Time Settings on page 5-4).

### CAUTION

---

Use DEXIS Sensor Holders and Positioning System only.

---

3. Assemble the rings, bars and holders so that the entire sensor area is visible through the ring. This applies to horizontal or vertical Bitewing holders, and Periapical (anterior or posterior) holders.
4. Apply a new disposable sheath to the image sensor, making sure that it covers the portion of the cable which may come into contact with the patient.

### CAUTION

---

When twisting the disposable covers around the sensor wire, be sure not to twist the wire.

---

5. Attach the sensor to the positioner. The sensor will only fit into the holder one way.
6. Position the sensor appropriately in the part of the mouth to be X-rayed (refer to specific instructions for use).

**Important:** The sensor must be positioned with the sensitive area facing the source of the radiation. The sensitive side of the sensor is marked with the DEXIS logo.

**NOTE:** Use of a sensor positioning device is recommended to ensure that the sensor is positioned at right angles to the radiation beam and parallel to the tooth. Positioning devices are supplied with the DEXIS system.

## Portability

DEXIS Sensor can be easily moved from one dental chair to another. Once disconnected from the USB port of the computer, the sensor may be positioned near a second dental chair and connected to a USB port present nearby.

The USB technology allows easy connection and disconnection of DEXIS Sensor even if the computer is switched on, unless the PC presents particular limitations. No activation or deactivation procedure is necessary with DEXIS Sensor, just plug in and out.

## Hygiene

### **CAUTION**

---

After each patient, the sheath must be removed and the sensor must be cleaned and disinfected by CaviWipes; holders should be disinfected in an autoclave or with a cleaning agent after each patient.

---

Store disposable sanitary sheaths in a clean, dry place not exposed to sunlight or UV rays.

Make sure that used sheaths are disposed of as infected waste which is potentially biologically hazardous.

**Important:** Disinfect the image sensor before its first use and whenever there is a risk of contamination.

Please follow the disinfecting and cleaning instructions in order to avoid damaging the sensors. Disinfection of the X-ray image sensor at the end of each day is recommended even when sheaths are used.

**WARNING**

Do not autoclave sensor. Autoclave sterilizers will permanently damage the DEXIS Sensor.

---

## Recommended Disinfecting Solutions

Decontaminate the sensor, cable and holders in accordance with CDC – or your country's (e.g., OSAP) standards – for infection control. DEXIS recommends wiping the DEXIS Sensor with one of the following cleaning agents.

- CaviCide® (Manufactured by Metrex)
- CaviWipes™ (Distributed by Kerr)
- ProSpray™ (manufactured by Certol)
- Sani-Cloth® Plus, HB (manufactured by PDI)

**DO NOT:**

- Clean the sensor using non appropriate tools.

## Maintenance

DEXIS Sensor does not require any special maintenance other than regular cleaning and disinfection. It is recommended use DEXIS Software Manager for automatic software updates.

Go to [dexis.com/downloads](http://dexis.com/downloads) for more information on using DEXIS Software Manager. DTX Go may also be used for managing your software updates.

## Image Quality Assurance

Image quality of the DEXIS Sensor depends on several factors:

- the quality of the X-ray source (kV, focal spot size, distance)
- the alignment of the X-ray source to the anatomic region
- the applied X-ray dose / exposure time
- the settings of the computer monitor

It is recommended that you establish a procedure for periodic review of the image quality. If image quality is not satisfactory, or degrading, please check the contributing system parts as outlined below:

## X-ray Dose / Exposure Time Settings

The DEXIS Sensor has been designed to be used with a wide range of dose settings between 20 uGy and 2,000 uGy, to allow adjustments of the dose to the specific diagnostic task, and to compensate for under- and over-exposure. As a general recommendation, start by using dose / exposure time settings recommended by your X-ray source manufacturer for digital X-ray sensors.

The DEXIS Sensor can be used with much lower dose / exposure time settings than radiographic film. Nevertheless, low dose operation of a digital sensor in general can result in a grainy image appearance. If the image of the DEXIS Sensor appears grainy, increase the dose settings. If you obtain good results with a particular setting, you might try a lower dose setting to see if you still gain good results.

The DEXIS Sensor can work at higher dose /exposure time settings, if needed for a particular diagnostic task. At high dose / exposure time settings, it might not be possible to distinguish air from soft tissue. Such regions may appear overexposed. If air and soft tissue regions in the image appear overexposed, reduce the dose settings.

Remember, as with standard film, you will need to adjust the duration setting to compensate for tooth type (central incisor to molar) and patient body type (larger adult to small child). Ultimately, the settings you choose should be what you consider to be suitable for your diagnostic needs.

## X-ray Sharpness / Contrast

Many parts of the X-ray imaging system contribute to the sharpness and contrast of the image. It is recommended that you use a dental phantom for a periodic (at least annual) assessment of the image quality by performing a side-by-side comparison of an initial dental phantom image to a current dental phantom image to ensure that the images have sufficient quality for diagnostic purposes in your application.

Dental phantom images should be acquired with fixed settings for X-ray dose (kV, mA, distance), and a fixed and reproducible alignment of the X-ray source and the dental phantom to the DEXIS Sensor.

**NOTE:** Contact DEXIS technical support for recommendations about commercially available dental phantoms.

## Display Image

Refer to the software manual for guidance on how to ensure good display settings and image display properties. It is recommended to view images on a factory calibrated LCD/LED monitor with at least 2,560 x 1,440 pixel resolution at 0.25 to 0.26 dot pitch or better and 32bit true color, in native resolution.

# Chapter 6 Specifications and Standards

## CAUTION

US Federal law restricts this device to sale by or on the order of a dentist or other licensed practitioner.

## Sensor Specifications

Sensor Dimensions	Titanium: 39.9 x 29.8 x 8.6 mm IXS Size 1: 37.0 x 25.2 x 8.3 mm IXS Size 2: 43.3 x 30.4 x 8.5 mm
Sensor Image Area	Indirect converting dental intraoral X-ray sensors Titanium: 1692 by 1324 pixels IXS Size 1: 1539 by 1026 pixels IXS Size 2: 1842 by 1324 pixels, <i>All: 19.5 μm pixel size</i>
Resolution	20+ visible lp/mm
X-ray Parameters	Sensor can be used with dental X-ray generators in the range of 60 to 70 kV; at minimal 40 μGy incident dose
Software Architecture	Operating Systems Supported: <ul style="list-style-type: none"> <li>• Microsoft Windows 7 with current service pack</li> <li>• Microsoft Windows 8.1 with current service pack</li> <li>• Microsoft Windows 10 with current service pack</li> </ul>
Electrical Rating	DC 5V, 0.5 W <span style="float: right;">— — — - - - -</span>
Connection to PC	USB 2.0 compliant / USB 3.x compatible
Protection Against Shock	Type BF Applied Part 
Mode of Operation	Continuous
Method of Sterilization	Sensor not suitable for sterilization

Environmental Conditions	Humidity	Air Pressure	Ambient Temperature
<b>USAGE</b> DEXIS Sensor is not suitable to be operated in oxygen rich and/or explosive environments	30% to 95%	700 to 1060 hPa	5° to 30°C
<b>TRANSPORTATION AND STORAGE</b> Transport in the supplied protective package	10% to 95%		-40° to 70°C
EU Classification	Class IIa medical device according to MDD 93/42/ECC		
Protection against water/matter – IP 68			

## Compatibility with Radiographic Generators

DEXIS Sensor is generally compatible with any dental X-ray unit and generator capable of supplying the required range of exposure times and doses.

Follow the instructions of the X-ray generator to set the desired dose range.

The radiation of X-ray tubes is controlled by the settings of:

- Exposure time (msec) or pulses
- Voltage (kV or kVp)
- Current (mA)

Some controls allow for the modification of all of the above, some have fixed settings for current and voltage.

### Distance of X-ray Source to Sensor

There is a correlation between the distance of the X-ray source to the sensor and the dose received by the DEXIS Sensor. The radiation that reaches the sensor decreases with the square of the distance. If you double this distance, you receive only 1/4 of the radiation.

# Chapter 7 Troubleshooting

This section supplies information on some simple tests which the user may perform in the event of malfunctioning. Refer to the PC manual and the software manual for information on other types of malfunctions.

## System Does Not Acquire X-ray Images

1. Check the DEXIS Sensor USB connection to the PC.
2. Check that the software drivers are installed correctly.
3. Make sure that the active side of the DEXIS Sensor is facing the source of the X-ray and the active area is aligned with the X-ray beam.
4. Check the X-ray exposure settings and ensure that the X-rays were emitted.
5. Contact Technical Support at 888-883-3947.

# Appendix A Accessories

	Part	Description	
1	1.013.3253	Universal Sensor Sheath (Qty. 100)	
2	1.013.3254	Universal Sensor Sheath (Qty. 1000)	
3	1.013.3252	Positioning System Sheath (Qty. 500)	
4	1.008.2665	USB 15' Extender Cable	
5	1.014.1145	USB 3' Extender Cable	
6	1.013.3841	Universal Sensor Garage Kit	
7	1.013.8980	USB Hub (StarTech #ST4300USB3)	

# Appendix B EMC Information

The DEXIS sensor is, like any electronic medical device, subject to electromagnetic interactions with other electronic devices. The information in this chapter addresses this issue.

The EMC information in this chapter is provided for the medical system established by connecting the DEXIS sensor to a computer. This computer must be compliant with IEC 60950-1 (if located outside the patient environment) or IEC 60601-1 (if located inside the patient environment). Please consult the documentation of the computer for complete EMC information.

The DEXIS sensor is suitable for use in hospitals except for near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbance is high.

**Important:** Portable/mobile radio frequency communications equipment can affect the function of the DEXIS sensor as well as any other electronic medical equipment. This affect may result in Image failure or image quality degradation.

DEXIS is a USB compliant device and shall be used with USB compliant cables suitable for high speed/USB 2.0 cables. Such cables are either marked "USB 2.0" or "USB high speed." USB certified hubs can be used to extend the distance to the USB host/computer. The length of the cable connection to the hub or between hubs shall not exceed 5 m.

DEXIS Sensor is a USB compliant device and shall be used with USB compliant cables suitable for high speed/USB 2.0 cables. Such cables are either marked "USB 2.0" or "USB high speed." USB certified hubs can be used to extend the distance to the USB host/computer. The length of the cable connection to the hub or between hubs shall not exceed 5 m.

## **WARNING**

---

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

---

**WARNING**

---

Using non-USB compliant cables or hubs, or exceeding the maximum count of USB hub devices for extending the distance, can degrade the immunity of the DEXIS sensor to electromagnetic fields or increase the emission of electromagnetic fields from DEXIS Sensor. For USB Hub, use 1.013.8980.

For cable extensions, use USB Cable Extender (3 ft) 1.014.1145 or the USB Cable Extender (15 ft) 1.008.2665.

---

**WARNING**

---

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

---

**WARNING**

---

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the DEXIS sensor, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

---

<b>Guidance and Manufacturer's Declaration - Electromagnetic Emissions</b>		
The DEXIS sensor, used with a compliant computer, is intended for use in the electromagnetic environment specified below. The customer or the user of the DEXIS sensor should assure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The sensor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The sensor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class B (*)	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies (*)	

(\*) Computer used with the DEXIS sensor must meet this rating.

<b>Guidance and Manufacturer's Declaration - Electromagnetic Immunity</b>			
The DEXIS sensor, used with a compliant computer, is intended for use in the electromagnetic environment specified below. The customer or the user of the DEXIS sensor should assure that it is used in such an environment.			
<b>Immunity Test</b>	<b>IEC 60601 Test Level</b>	<b>Compliance Level</b>	<b>Electromagnetic Environment - Guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	Complies	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines	Complies (*)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1$ kV differential mode $\pm 2$ kV common mode	Complies (*)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% $U_T$ , 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% $U_T$ ; 1 cycle 70% $U_T$ ; 25/30 cycles for 50 Hz and 60 Hz respectively Single phase: at 0° 0% $U_T$ , 250/300 cycle for 50 Hz and 60 Hz respectively Single phase: at 0° <b>VDI specs per 3<sup>rd</sup> edition:</b> <5% $U_T$ (95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycle 70% $U_T$ (30% dip in $U_T$ ) for 25 cycle <5% $U_T$ (95% dip in $U_T$ ) for 5 s	Complies (*)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the DEXIS sensor requires continued operation during power mains interruptions, it is recommended that the DEXIS sensor be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50 Hz or 60 Hz	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE 1 $U_T$ is the a.c. mains voltage prior to application of the test level.			

(\*) Computer used with the DEXIS sensor must meet this rating.

### Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The DEXIS sensor, used with a compliant computer, is intended for use in the electromagnetic environment specified below. The customer or the user of the DEXIS sensor should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
<p>Conducted RF IEC 61000-4-6</p>	<p>AC Mains 3V, 0.15 MHz - 80 MHz 6 V in ISM band between 0.15 MHz and 80 MHz 80% AM at 1 kHz</p> <p>SIP/SOPS: 3V, 0.15 MHz - 80 MHz 6 V in ISM band between 0.15 MHz and 80 MHz 80% AM at 1 kHz</p>	<p>Complies (*)</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the DEXIS sensor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p><math>d = 3.5 / (V1) \sqrt{P}</math>    150 kHz to 80 MHz</p> <p><math>d = 3.5 / (E1) \sqrt{P}</math>    80 MHz to 800 MHz</p> <p><math>d = 7 / (E1) \sqrt{P}</math>    800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p><math>V_1 = 3 \text{ Vrms}</math> <math>E_1 = 3 \text{ V/m}</math></p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>
<p>Radiated RF IEC 61000-4-3 ed3.0 (with A1:2007 +A2:2010)</p>	<p>3 V/m 80 MHz to 2,7 GHz 80% AM at 1 kHz</p>	<p>Complies (*)</p>	<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  </div>

NOTE 1 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DEXIS sensor is used exceeds the applicable RF compliance level above, the DEXIS sensor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the DEXIS sensor.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

<b>Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the DEXIS Sensor</b>			
<p>The DEXIS sensor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the DEXIS sensor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the DEXIS sensor as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter in Watts	Separation distance according to frequency of transmitter in meters		
	150 kHz to 80 MHz $d = 3.5 / (f/1) \sqrt{P}$	80 MHz to 800 MHz $d = 3.5 / (f/1) \sqrt{P}$	800 MHz to 2,5 GHz $d = 7 / (f/1) \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance <math>d</math> in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

# Appendix C Product Symbols

Symbol	Title of the Symbol	Reference Number	Standard Containing the Symbol	Function / Description per Standard	Manufacturer Interpretation
	Type BF applied part	5333	IEC 60417	To identify a type BF applied part complying with IEC 60601-1.	To identify a type BF applied part complying with IEC 60601-1.
	Direct Current	5031	IEC 60417	To indicate on the rating plate that the equipment is suitable for direct current only; to identify relevant terminals.	Direct Current
	Intertek ETL Listed (Canada & USA)	N/A	N/A	The Marks are evidence that the model/product conforms to applicable standards and that there is a program of ongoing factory inspections.	This ETL listed mark guarantees that Intertek has certified the product described herein under control number 5009700 to be in compliance with the applicable regulations. Intertek is: <ul style="list-style-type: none"> <li>• a Nationally Recognized Testing Laboratory by the Occupational Safety and Health administration (OSHA) in the United States.</li> <li>• a Certifying Body in Canada by the Standards Council of Canada.</li> </ul>
	CE Mark	N/A	N/A	The product meets all the legal requirements for CE marking and can be sold throughout the European Economic Area.	The CE symbol ensures that the product specified herein meets the provisions of European Council Directive 93/42 EEC concerning medical devices.
	General Warning Sign	W001	ISO 7010	To signify a general warning	Warning
	Warning; Radioactive material or ionizing radiation	W003	ISO 7010	To warn of radioactive materials or ionizing radiation	Warning
	Caution	0434B	ISO 7000	To indicate that caution is necessary when operating the device or control close to where the symbol is placed, or to indicate that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.	Please refer to the written instructions of this manual.

Symbol	Title of the Symbol	Reference Number	Standard Containing the Symbol	Function / Description per Standard	Manufacturer Interpretation
	Do not reuse	1051	ISO 7000	To indicate that the item is for single use only and must not be used more than once, for example on packages of medical disposables.	Indicates the product should be used only once.
	Marking of electrical and electronic equipment	N/A	BS EN 50419	This symbol on the products and/or accompanying documents means that used electrical and electronic products should not be mixed with general household waste.	Please contact your dealer or supplier for further information about product disposal at the end of the product's lifetime.
	Manufacturer	3082	ISO 7000	To identify the manufacturer of a product. This symbol shall be used filled in all applications to differentiate it from ISO 7000-2497.	Manufacturer
	Date of manufacture	2497	ISO 7000	To indicate the date on which a product was manufactured.	Date of manufacture
	Catalogue number	2493	ISO 7000	To identify the manufacturer's catalogue number, for example on a medical device or the corresponding packaging. The catalogue number shall be placed adjacent to the symbol.	Catalogue number
	Serial number	2498	ISP 7000	To identify the manufacturer's serial number, for example on a medical device or its packaging. The serial number shall be placed adjacent to the symbol.	Serial number
	Operator's manual; operating instructions	M002	ISO 7010	To identify the location where the operator's manual is stored or to identify information that relates to the operating instructions. To indicate that the operating instructions should be considered when operating the device or control close to where the symbol is placed.	Follow Operating Instructions for Use
<b>IP68</b>	N/A	N/A	N/A	N/A	Protection against dust and continuous immersion in water.

Symbol	Title of the Symbol	Reference Number	Standard Containing the Symbol	Function / Description per Standard	Manufacturer Interpretation
	Atmospheric pressure limitation	2621	ISO 7000	To indicate the acceptable upper and lower limits of atmospheric pressure for transport and storage.	Atmospheric pressure limitation
	Humidity limitation	2620	ISO 7000	To indicate the acceptable upper and lower limits of relative humidity for transport and storage.	Humidity limitation
	Temperature limit	0632	ISO 7000	To indicate the maximum and minimum temperature limits at which the item shall be stored, transported or used.	Temperature limitation
	This way up	0623	ISO 7000	To indicate correct upright position of the transport package.	This end up
	N/A	N/A	ISO 7000	N/A	Stacking unit
	Keep away from rain	0626	ISO 7000	To indicate that the transport package shall be kept away from rain and in dry conditions.	Keep dry
	Fragile; handle with care	0621	ISO 7000	To indicate that the contents of the transport package are fragile and the package shall be handled with care.	Fragile handle with care
<b>R<sub>x</sub> Only</b>	Prescription Device	801.109	FDA Title 21	Caution: U.S. Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner	Clinical Use Only
	The name and address of the authorized representative in the European community.				

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

1. Connect your sensor to your computer's USB port to start the automatic download of the calibration files.

NOTE: Right click this icon to open the *Intra-Oral Sensor Dashboard*.

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**NOTA:** DO NOT unplug the sensor while the calibration files are being downloaded.

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**Answer:** Yes. Not doing so may cause the software and drivers to not operate correctly. The following folders need to be added for exceptions/exclusions depending on which drivers and imaging software you are installing.

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**PaloDEx Group Oy**