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SONOREX DIGITEC und SONOREX DIGITEC-RC

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Operating Instructions

SONOREX DIGITEC und SONOREX DIGITEC-RC

ultrasonic cleaners for aqueous cleaning agents

Cleaning units, accessories and cleaning agents have to be strictly used according to the instruction manual or the corresponding product information.

Keep cleaning units and cleaning agents away from children and from persons not instructed for correct use by means of these instructions.

In case of improper use, the manufacturer does not take any responsibility for safety and proper functioning.

Unauthorized changes/modifications will result in expiration of the CE-conformity.

If service is required, please contact your dealer or the address mentioned below.

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Schutzklasse: Klasse I

Contents

| | | |
|-------|---|-----------|
| 1 | Strictly observe the following general instructions for safety and correct operation! | 18 |
| 2 | Start | 19 |
| 2.1 | Degassing | 19 |
| 2.2 | Cleaning | 19 |
| 2.3 | Explanation of operating elements and start-up | 21 |
| 2.4 | Operating modes | 22 |
| 2.4.1 | Operation with heating (units with heating) | 22 |
| 2.4.2 | Ultrasound operation without heating | 22 |
| 2.4.3 | Ultrasound operation with heating | 23 |
| 2.5 | Leads for operation and special funktions | 24 |
| 2.6 | Test routine for operating functions | 25 |
| 2.7 | Emptying the ultrasonic tank | 25 |
| 3 | Additional recommendations | 26 |
| 3.1 | Cleaning agents | 26 |
| 3.2 | Cleaning and care | 26 |
| 3.3 | Treatment of medical instruments | 27 |
| 4 | Disfunctions | 27 |
| 5 | Dosing table for SONOREX – compact ultrasonic cleaners | 28 |
| 6. | Technical data | 29 und 44 |

1 Strictly observe the following general instructions for safety and correct operation!

- Keep cleaning units and cleaning agents away from children and also away from persons not instructed for correct use by means of these instructions.
- Do not immerse parts of the body (hand, feet) nor creatures or plants into the ultrasonic tank. During ultrasonic cleaning do not dip your hands into the cleaning liquid.
- Avoid any contact of your eyes and/or your skin with the cleaning liquid.
- Do not ingest nor breathe in the cleaning liquids.
- Never run the unit without liquid!
- Do not fill hot water into the ultrasonic tank. The filling temperature must not exceed 50 °C.
- Do not put parts to be cleaned into the tank with aggressive soiling such as acids, chloride ions or similar substances!
- Do not fill the stainless steel oscillating tanks with aggressive cleaning agents such as acids or chemicals which contain or release chloride ions (some disinfectants, dish-washing detergents, household cleaners, saline solutions).
- Do not fill the stainless steel oscillating tanks with combustible liquids.
- Use aggressive cleaning agents only in inset beakers or non-perforated plastic insert tanks.
- After removing the packing, check the unit's casing and control elements for any signs of damage. In such case, do not connect the unit to the electrical power supply, but report the damage immediately and in writing to the carrier and your supplier. Save the original packing.
- Place the unit in a horizontal position on a dry, solid base.
- Plug the unit only to a socket with grounded safety contact.
- Keep the unit's surface and control elements clean and dry!
- Replace defective parts only with original **SONOREX** parts!
- Ultrasonic cleaning produces cavitation noise. If working within distances of 2 meters from the unit for longer periods of time, be sure to wear ear plugs to avoid damage to your hearing!
- Do not place anything on the bottom of the oscillating tank. Use appropriate **SONOREX** accessories (insert basket, positioning lid with insert vessel, etc.).
- The ultrasonic energy in the oscillating tank heats up the cleaning liquid without additional heating! Check the temperature regularly when cleaning temperature-sensitive parts.
- Do not move or carry units filled with liquid.

Principle of ultrasonic cleaning

PZT oscillators underneath the tank convert electrical energy into mechanical vibrations.

SONOREX units produce oscillations of 35 kHz in the cleaning liquid. Tiniest vacuum bubbles are formed which then implode (cavitation). This cavitation principle removes dirt and grease from the deepest pores - "electronic scrubbing" - at all points that come into contact with the cleaning liquid.

2 Start

- Make sure that the ball valve of the unit is closed!
- Fill the oscillating tank up to the filling mark.
- Connect the unit to a grounded socket.

2.1 Degassing

Degassing the cleaning liquid enhances the cleaning effect - see also lead for DEGAS function under chapter 2.5.


Cleaning and contact liquids contain dissolved gases (for example, oxygen). Sonicate cleaning and contact liquid, which have just been poured into the oscillating tank or have been left in it for a longer period of time, for approx. 5 to 15 minutes before starting the cleaning process. Degas oscillating tanks with a volume of more than 10 l for roughly 30 minutes. During the degassing process the cavitation noise changes. As the loud degassing noise decreases towards the end of the process, the unit appears to run more quietly. Lower noise levels indicate the end of the degassing process and an improved cleaning impact. They do not mean a reduction of ultrasonic energy.

2.2 Cleaning

- Be sure that the objects in the tank are completely covered with liquid.
- Normally the objects may be cleaned **directly** in the ultrasonic tank with an insert basket.
- In the following cases clean objects **indirectly** in inset beakers or non-perforated insert tanks to protect the stainless steel oscillating tank:
 - when chemically aggressive cleaning liquids (for example acids, etc.) are used.
 - when chemically aggressive substances (for example when cleaning racks from film developing machines) or
 - when abrasive substances (for example polishing paste, quartz, sand) are to be removed.


DIRECTLY in the oscillating tank:

- Fill the oscillating tank to 2/3 with cleaning liquid (TICKOPUR, STAMMOPUR see dosage on label). Improperly filled units can suffer damage (see point 5 – dosing table).
- Use always an insert basket or utensil holder.
- Parts must not be placed on the tank bottom (Exemption basket K6 and SH 7).


INDIRECTLY in inset beakers or non-perforated insert plastic tub
⇒ **SONOREX** Accessories

- Fill the oscillating tank with a tenside-containing contact liquid (for example TICKOPUR R 33 or STAMMOPUR RD 5 diluted with water) to achieve a good transmission of ultrasound.
- Fill the inset beakers or insert plastic tub with cleaning liquid (for example TICKOPUR, STAMMOPUR).
- Suspend inset beakers in the positioning cover, do not place them on the tank bottom.
- Inset beakers or insert plastic tub must be immersed in the contact liquid to a depth of at least 2 cm. See to the tank filling to 2/3.
- Several inset beakers with different cleaning liquids can be sonicated at the same time.


When inserting the objects to be cleaned, be sure to observe the following:

- Place the dirtier side towards the bottom.
- Do not stack - ultrasonic energy will be absorbed.
- Sensitive pieces may not be allowed to touch one another. Avoid damage by using plastic insets, plastic insert baskets or the knobbed silicone mats in insert basket.
- Be sure that there are no air bubbles in hollow spaces (for example pocket holes).
- Check the temperature when cleaning heat-sensitive objects.

2.3 Explanation of operating elements and start-up

example: unit with heating

| symbol | operation |
|--------|--|
| | Switch the unit on or off (standby). After switching on the unit, the heating is always off. Time scale shows last value. Actual temperature is indicated through blinking LED. |
| | Setting of target temperature through keypress. Permanent pressing (> 2 sec) switches off the heating. |
| | Time setting through keypress. Permanent operation = ∞. |
| | Start and Stop of ultrasound through keypress. |

2.4 Operating modes

2.4.1 Operation with heating (units with heating)

- Switch on the unit (press key On/Off) - green LED on the temperature scale lights up.
- Select target temperature – increase value gradually by pressing the Temperature key
target temperature = lighting yellow LED
actual temperature = blinking yellow LED
- When target temperature is higher than actual temperature, heating phase begins immediately. As long as heating is switched on, ultrasound is being switched on for 3 sec per min. See under 2.5 temperature – avoid retardation of boiling.
- As soon as the target temperature is reached, heating switches off = lighting LED.
- Switch off heating manually:
Press Temperature key longer than 2 sec – green LED lights or switch off the unit with On/Off-key (all LEDs off).
Last selected temperature will not be stored.

Lead: Target temperature can be modified at any time independently from ultrasound operation.

2.4.2 Ultrasound operation without heating

- Switch on the unit (press key On/Off) – yellow LED of last selected time on the time scale lights up.

Start without modification of time:

- By pressing the key Start/Stop, ultrasound is being switched on.
This is indicated on the time scale by means of a running light of the balance time.
- By pressing the key Time, ultrasound output can be interrupted (pause).
- By pressing the key Time again, ultrasound output will be continued (balance time).

Start with modification of time:

- Select time - by pressing the key Time, the duration will be gradually increased.
Numeral (1 ... 30) = lighting, yellow LED
Continuous operation (∞) = lighting, green LED
- By pressing the key Start/Stop, ultrasound is being switched on.
This is indicated on the time scale by means of a running light of the balance time.
- By pressing the key Time, ultrasound output can be interrupted (pause).
- By pressing the key Time again, ultrasound output will be continued (balance time).
- After run off of selected time, the unit stops and LED of previously selected time lights up.
- Switch off the unit by pressing the key On/Off.
The time last selected will be stored.

Attention! Observe leads for operation and special functions – chapter 2.5

Cleaning time

Cleaning time must be kept as short as possible. Depending on the degree of soiling and the cleaning liquid used (see the data on its label), the cleaning time varies between 1 and 15 minutes.

2.4.3 Ultrasound operation with heating

- Switch on the unit (press key On/Off) – green LED on the temperature scale and yellow LED of the last selected time on the time scale lights up.
- Select target temperature – increase value gradually by pressing the Temperature key
target temperature = lighting, yellow LED
actual temperature = blinking, yellow LED
- When target temperature is higher than actual temperature, heating phase starts immediately. As long as heating is switched on, ultrasound is being switched on for 3 sec per min. See under 2.5 temperature – avoid retardation of boiling.
- As soon as the target temperature is reached, heating switches off = lighting LED.
- Switch off heating:
Press Temperature key longer than 2 sec – green LED lights up.

Start without modification of time:

- By pressing the key Start/Stop, ultrasound is being switched on.
This is indicated on the time scale by means of a running light of the balance time.
- By pressing the key Time, ultrasound output can be interrupted (pause).
It is possible to modify the target temperature by pressing the key Temperature.
- By pressing the key Time again, ultrasound output will be continued (balance time).

Start with modification of time:

- Select time - by pressing the key Time the duration will be gradually increased.
Numeral (1 ... 30) = lighting, yellow LED
Continuous operation (∞) = lighting, green LED
- By pressing the key Start/Stop, ultrasound is being switched on
This is indicated on the time scale by means of a running light of the balance time.
- By pressing the key Time, ultrasound output can be interrupted (pause).
It is possible to modify the target temperature by pressing the key Temperature.
- By pressing again the key Time, ultrasound output will be continued (balance time).
- After run off of selected time, the unit stops and LED of previously selected time lights up.
- Switch off the unit by pressing the key On/Off.
The last selected time will be stored.

Attention! Observe leads for operation and special functions – chapter 2.5

2.5 Leads for operation and special functions

Important leads regarding temperature

- Warmed-up liquids improve the cleaning efficiency.
- However, the effect of ultrasonic cavitation decreases at higher temperatures¹. Experience has shown that best cleaning results are achieved by using an aqueous solution with a bath temperature of 50 to 60 °C.
- Ultrasonic energy also heats the cleaning agent. Continuous irradiation and/or the covering of the oscillating tank may cause the temperature of the liquid to rise more quickly and exceed the default setting.
- **Do not** heat disinfectants.
- When reaching a critical temperature, the ultrasonic power is reduced to protect the electronic components inside the ultrasonic device. This way, the internal temperature is prevented from further rising.

Temperature scale:

- LED of target temperature lights up permanently.
- LED of actual temperature is blinking. Above the adjusted target temperature LED is blinking slowly (= heating is off), below the target temperature LED is blinking quickly (= heating is on). When the actual temperature is equal to the target temperature, only the LED of the target temperature lights up. Display accuracy: $\pm 2,5$ °C.
- The LED „!“ will blink red when temperature exceeds 80 °C.

Temperature - avoid retardation of boiling:

- After switching on the unit, the temperature pre-selection is basically adjusted to „0“ (off).
- When a temperature is set, the unit will immediately try to reach the selected temperature. When heating (without ultrasound) switched on ultrasound will be automatically switched on once during 3 sec every minute for mixing the liquid and thus avoiding retardation of boiling.
 - This function is always active with temperatures > 60 °C and cannot be deactivated!
 - The function can be activated for temperatures < 60 °C. It must be activated again after each switch-on: Activation of the function:
Press the Temperature key while at the same time switching on the unit (On/Off).

Time scale / duration:

- LED of selected time is lighting constantly, in case of continuous operation (∞) the top (green) LED lights up.
- Ultrasound activity is indicated through LED blinking in second cycle (running light downwards). Running light always starts at the LED for remaining time.

Time scale / DEGAS:

- For degassing (DEGAS) before cleaning, switch on the DEGAS function. If necessary, select required duration with the key Time. Then press the key Start/Stop for at least 2 sec. Early switch-off through re-pressing the Start/Stop-key.
During degassing, the top green LED is blinking in addition to the time lapse (▬▬▬)
- Switch between ultrasound and degas: if the Start/Stop-key is pressed a long time with ultrasound on, the ultrasound is first turned off and then reactivated with the degas function after about 2 sec.

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Time scale / continuous operation:

- In order to avoid inadvertent switching on of continuous operation, continuous operation can be deactivated:
 - Pull the mains plug.
 - Press the Time key and simultaneously pocket the mains plug.
The yellow LED 1 min lights up for confirmation.The function is again activated in the same way. The green LED continuous operation (∞) lights up for confirmation.

Attention:

For safety reasons the unit will be automatically switched off, if no key has been pressed longer than 12 hours.

2.6 Test routine for operating functions

A test routine can be started in order to carry out an internal function test:

For this, the unit must be switched off. With the key Start/Stop pressed switch on the unit with the key On/Off.

One after the other, all LEDs will light up for one third second. Afterwards, the last selected values will reappear.

With this, the test has been completed successfully.

In case of deviations, the unit has to be sent in for inspection/repair.


2.7 Emptying the ultrasonic tank

Pull the mains plug.

Do not place the unit into a sink.


Units with drain

- Ball valve
Turn the lever to the outflow direction.

Units without drain

Pour out the liquid over the rear left tank corner. Avoid splashing from below. Pay attention to the cable connection.

3 Recommendations

3.1 Cleaning agents

- Only use approved and especially developed cleaning agents, e.g. TICKOPUR, TICKOMED, STAMMOPUR in order to obtain a high cleaning efficiency. Tap water without any additives does not clean sufficiently.
- When using the cleaning agents you must always strictly observe the safety instructions in the product information.
- Keep cleaning agents away from children and from persons not instructed for correct use by means of the product information.
- Do neither ingest nor inhale the cleaning agents. Avoid any contact of your eyes and/or skin with the cleaning agents.
- Replace used liquids and disinfectants, do not “freshen them up” with an added dose.
- Layers of dirt on the tank bottom or in the insets (abrasives and polishing agents, for example, as well as other deposits) reduce the cleaning efficiency. Empty the oscillating tank and clean it regularly. Replace the cleaning liquid.
- The cleaning and contact liquids in the ultrasonic oscillating tank may not exceed a maximum operating temperature of 100 °C.
- Distilled or deionized water without additives is unsuitable for ultrasonic cleaning directly in the oscillating tank. They lack the necessary cavitation properties to form a uniform ultrasonic field. Irregular cavitation can damage the oscillating tank or the objects being cleaned. (You may use distilled or deionized water in the inset beakers and non-perforated insert tubs.)
- Be careful when working with aggressive cleaning agents in insets or non-perforated insert tanks; avoid splashing into the contact liquid or onto stainless steel surfaces; if necessary, replace the contact liquid immediately, clean the surfaces and wipe them dry.

The manufacturer does not take any responsibility for damage of equipment or damages at the parts to be cleaned caused by inappropriate cleaning chemicals.

3.2 Cleaning and care

- Be sure to pull the plug before cleaning the unit.
- Do not shower the units or immerse them in water; instead, wipe them with a damp cloth, in- and outside.
- Rinse out the stainless steel tank frequently and wipe it dry.
- Do not use steel wool, scrapers or graters.
- If rings should appear in the stainless steel tank after a longer period of operation, remove them with a customary stainless steel cleaner without scouring additives.
- Metal parts remaining on the stainless steel surface, as well as particles of rust, can penetrate the passive protective layer of the stainless steel, thereby “activating” it and causing it to rust. This rust provokes pittings on the stainless steel surface. Therefore, remove all remaining metal parts such as screws, metal chips and the like; remove minor rust spots immediately with a soft cloth and customary stainless steel cleaners without scouring additives.

SONOREX units do not require maintenance.

3.3 Treatment of medical instruments

- The ultrasonic cleaning unit is classified as medical device, class II
- When disinfecting and cleaning contaminated medical instruments in an ultrasonic bath, the hygienic safety is important after the treatment. In case of improper and irregular disinfection and cleaning a micro-biological contamination is possible caused by settlement of microorganism, especially at the tank rim or in the drain section. Cross infection might be a result. Therefore, the tank and housing surfaces must be regularly disinfected and cleaned. According to the hygienic plan the user must carry out regular disinfection and cleaning with a surface disinfecting agent certified by VAH, or another efficiently classified agent.
- The instruments may not be placed on the tank bottom. The basket avoids damages at the instruments and at the tank bottom. Do not staple instruments, overloading reduces the cleaning efficiency. Pliers and scissors must be opened completely or detached, if necessary. Instruments must be completely covered with cleaning liquid, air has to escape from cavities and hoses.

Remark: Supplementary application information has to be respected.

4 Disfunctions

- Oscillations are too weak and non-uniform, the oscillation noise is too loud; the cleaning effect is insufficient:
 - Has the liquid been properly degassed? ⇒ sonicate for 15 minutes.
 - Is the unit overloaded? ⇒ remove some parts.
 - Non-uniform noise (wobbling) is not a defect ⇒ change the fluid level.
- Heating defect:
No objection to operate the unit without the heating.
- Minor erosion marks at the tank bottom? ⇒ common wear and tear effect.
The unit is faultless.
- The mains fuse can be released. ⇒ The unit is equipped with an overvoltage arrestor which protects it from voltage peaks; the mains fuse can be released during overcharge.
- Repairs should only be carried out by authorised service.
- Replace defective parts only by original **SONOREX** parts.

Please inform us about disfunctions in writing.

5 Dosing table for SONOREX – compact ultrasonic cleaners

You can order this dosing table (DIN A 4-size) free of charge or download it as PDF-file from the following internet address:

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Values have been rounded up mathematically.

number printed in normal type: water

number printed in **bold** type: concentrate

Fill the oscillating tank to 1/3 with water and add the appropriate quantity of concentrate. Then add water up to the filling mark and sonicate for approx. 15 minutes for degassing.

or

According to the dosing table, fill the oscillating tank with water, add the appropriate quantity of concentrate, mix the solution for a short while and sonicate for approx. 15 minutes for degassing.

Calculation of other volumes

Example:

- 10 litres solution ready for use
- 2,5 % dosage of the concentrate

$$\frac{10 \text{ l} \times 2,5 \%}{100 \%} = 0,25 \text{ l concentrate}$$

$$10 \text{ l} - 0,25 \text{ l} = 9,75 \text{ l water}$$

6 Technical data

SONOREX compact units are RFI-proofed and C € marked.

Compliance with the limits according to EN 55011, EN 61000-6-1/08.2002 and EN 60601-1-2.

Nominal voltage: 230V~ 50/60 Hz, (115 V on request)
cable length 2 m
High frequency: 35 kHz
Degree of protection: IP 33, **SONOREX DIGITEC-RC** IP 23
Oscillating tank: stainless steel

Environmental conditions in accordance with EN 61 010-1 (IEC 1010-1)

Degree of soiling: 2 according to IEC 60664-1:11.02
Excess voltage category: II
Permissible ambient temperature: 5 - 40 °C
Permissible relative humidity up to 31 °C: 80 %
Permissible relative humidity up to 40 °C: 50 %
Dew not permitted.

Information for the use as a medical product

Denotation: ultrasonic cleaning device
UMDNS-nomenclature (ECRI): 14-263
Purpose determination: ultrasonic-disinfection and cleaning of medical and dental instruments
Classification (Medical Devices Directive 93/42/EEC, appendix IX): class I; active, not invasive, not implantable medical device
Type, model, serial number, year of manufacture: see type designation plate on rear of device

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