Anwendertipps User tips

Inhalt

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GUTTA MASTER®

There are numerous scientific studies supporting the carrier-based technology which is today proven and established. The GuttaMaster® system is a warm filling method which can be learned quickly and carried out easily.

Before using GuttaMaster® for the first time

Make yourself familiar with the GuttaMaster® filling process!

Before the first clinical use of the GuttaMaster® filling process, we recommend practising on plastic blocks (included in the Starter Kit) and extracted teeth, which can be prepared with a NiTi system such as Mtwo®, to make yourself familiar with this filling technology.

Phase 1

Prepare the root canal for the obturation with the GuttaMaster® system


It is important that you prepare the root canal to a conicity of 5 % or higher in the apical region. Even tight or curved root canals must be prepared to this standard and at least to the smallest verifier size. Thanks to the high flexibility of the Mtwo® system, a conicity of 6 % can already be obtained with a tip diameter of 20.

Next to an efficient preparation of the root canal, a thorough cleaning of the canal is an essential prerequisite for a successful obturation.

Tip: Use of optical magnifying aids is recommended, e.g. surgical loupe or dental microscope.

In order to select the right GuttaMaster® obturator size, we recommend using a verifier. Adjust obturation reference length (approx. 1 mm to working length) with the help of the ring markings and/or the stopper, for example by using the Minifix measuring gauge.

Adjust obturation reference length

In finding the right GuttaMaster® obturator, the goal is that the verifier exhibits a gripping effect approx. 1 mm to working length (obturation reference length) when extracted.

Procedure: In order to choose the right GuttaMaster® obturator size, insert the verifier into the root canal up to the obturation reference length. Start with a verifier corresponding to the tip diameter of the largest preparation instrument used at working length.

The following cases can occur:

A) The verifier goes to working length and exhibits a gripping effect. You have selected the right size. Kindly select the colour-corresponding obturator.

B) The verifier goes to working length but exhibits no gripping effect. Please select the next largest verifier until a gripping effect is achieved at obturation reference length. If you have to go two or more sizes higher, please reconsider the earlier preparation. If necessary, the root canal has to be extended to higher apical diameters.

C) The verifier clamps before obturation reference length. Kindly select the next smaller verifier until a gripping effect is achieved at obturation reference length.

Please note: In that case please reconsider the carried out preparation critically and verify that your preparation has extended to full working length.

Select the appropriate GuttaMaster® obturator

Note:
In what follows, “approx. 1 mm to working length” is described as the obturation reference length, as it is the decisive factor for determining the suitable obturator.

Note:
Depending on the situation, a larger or smaller verifier is selected until a verifier reaches the obturation reference length (approx. 1 mm to working length) and there exhibits a gripping effect. Please use the GuttaMaster® obturator in the corresponding colour for filling.

Phase 4

Disinfect GuttaMaster® obturator

Disinfect the GuttaMaster® obturator for at least one minute in a sodium hypochlorite or a high proof alcohol solution and then dry it.

Phase 5

Heat up the GuttaMaster® obturator in the GuttaMaster® oven

Note:

Before hanging up the GuttaMaster® obturator, the appropriate obturation reference length (approx. 1 mm to working length) should be set correctly.


Hang the GuttaMaster® obturator in the pulled out obturator holder of the GuttaMaster® oven.

Please ensure:

The silicon stop must be **underneath** the holder.

Select the appropriate size setting on your GuttaMaster® oven for the obturator to be heated. Press the  -key in front of the „LEFT“ or „RIGHT“ holder used.

If several root canals are to be filled, you can save time by using the left and right holder alternately (see GuttaMaster® oven Operating Manual, point 9).

Phase 6

Dry root canal, mix sealer and apply

Mix the root canal sealer (e.g. 2Seal). Dry the canal before carefully applying the sealer with a paper point according to the size and shape of your root canal preparation. With another paper point, apply a thin sealer layer to the root canal.

Then use a new paper point to remove excess sealer from the root canal. This procedure reduces the risk of possible overpressing.

Please ensure:

Kindly use a sealer suitable for warm obturation methods, e.g. 2Seal.


Phase 7

Fill root canal

Carefully take the GuttaMaster® obturator out of the GuttaMaster® oven holder. Take care that the gutta-percha is not inadvertently scraped off the plastic core, and insert the obturator with a slow, steady, apically directed movement up to the obturation reference length into the root canal. To prevent the carrier from bending, it is important that your movement is always directed along the longitudinal axis of the carrier.

On insertion, the gutta-percha flows apically in front of the carrier. This allows for a filling up to the apex, as well as a closing of side canals and isthmi. After having reached the obturation reference length, the pressure should be maintained for a few seconds, for the gutta-percha to cool off under light pressure (thus counteracting a temperature-related shrinking effect of the gutta-percha).

Phase 8

Separate core and handle

To conclude the root canal filling, the handle must be separated from the plastic core, for example by using a GuttaMaster® PREP® Bur (turbine 150,000-200,000 rpm without water cooling).

Phase 9

Remove excess gutta-percha

Note:

With a multiple root tooth, repeat the steps described above for all root canals.

Begin the filling of several canals in a tooth with the one which is most difficult to reach.

To avoid the unintended transfer of gutta-percha or sealer, it is helpful to initially seal off canals that will be filled at a later time with a foam pellet.

If necessary, compact the gutta-percha remaining around the core post, for example with a Machtou Plugger. If excess gutta-percha blocks the access to other canals, remove it with the help of a warm endodontic excavator.

After filling all root canals, sealer residue can be removed with a small sponge and alcohol.

First remove the coronal section of the filling with a GuttaMaster® PREPI® Bur, using a turbine or a contra-angle at a speed of 150,000-200,000 rpm. Do not spray and work with light pressure.

Then use the appropriate post drill to extend the canal to the size of the chosen post (e.g. DT Light® SL). You can then continue with your normal post setting routine.


Core points of the GuttaMaster® obturators have a longitudinal groove specially designed to ensure easy removal.

To remove the GuttaMaster® obturator, insert an appropriate rotating nickel-titanium instrument (e.g. the Mtwo® Retreatment file) in order to remove the gutta-percha and the plastic core. The gutta-percha can be softened up with a solvent (such as eucalyptol), thus facilitating the penetration of the revision instrument.

Alternatively, GuttaMaster® obturators can also be removed with a manual instrument such as a Hedstroem file.

If a root post
should be placed

If a post
should be revised

FAQ – Frequently Asked Questions

1. What is the GuttaMaster® system?

The GuttaMaster® system is based on the tried and tested carrier-based technology. A plastic carrier is coated in alpha gutta-percha with a low melting point. The GuttaMaster® system allows you to achieve a tight filling, including lateral canals and isthmi, efficiently and in just a few steps. Obturators and verifiers are adapted for preparation with the Mtwo® NiTi system.

2. Is the plastic carrier/verifier biocompatible?

Yes. The carrier/verifier material consists of a non-toxic polymer compound. This synthetic material is also used in implantology and is physically very well accepted.

3. Are the obturators sterilised?

As a rule, gutta-percha is not sterilised. Therefore the obturator should be disinfected before use. Disinfect the GuttaMaster® obturator for at least one minute in a sodium hypochlorite solution or in high proof alcohol and then dry it.

4. How important is it to dry the root canal before obturation?

This step is very important, as moisture in the root canal can cause an insufficient binding of sealer and gutta-percha resulting in a not tight enough filling.

5. What can I do if the stopper of the GuttaMaster® obturator is situated at the exact height of the oven holder?

Carefully remove the stopper, e.g. by using a scalpel and mark the obturation reference length, for instance with the help of a permanent marker.

6. Is a sealer needed for GuttaMaster®?

As with any root filling process that is based on gutta-percha as the core material, GuttaMaster® also requires a sealer for sealing purposes. In addition, the sealer improves the lubricity of the plasticised gutta-percha. The sealer smoothes out the unevenness and improves the connection between the gutta-percha filling and the dentinal wall.

Furthermore, using a sealer seems to reduce the temperatures acting on the outer radicular walls. The temperature-lowering effect of the sealer is highly appreciated, even if we know from several studies that the temperatures used with this filling method are not harmful to the periradicular structures.

7. Which type of sealer is recommended?

A sealer suitable for warm obturation methods should be used. We recommend the use of 2Seal. This sealer has been proven suitable for warm filling methods and furthermore is recommended as a filling material in the guidelines of the DGZMK (German Dental Association).

8. What amount of sealer is recommended?

Apply a small amount of sealer, preferably with a dry paper point around the root canal entrance, and spread the sealer with a stirring movement up to the apex. Please note that a thin sealer layer is completely sufficient. The root canal wall only has to be lightly coated.

9. Do I need to change my root canal preparation routine or my instruments when working with GuttaMaster®?

The GuttaMaster® system is tailored to preparation with the M_{two}® NiTi system. Therefore you should perform your preparation with M_{two}® and, as described in the user tips, prepare the canal form according to the GuttaMaster® filling method.

10. Why is a conicity of 5% or higher a prerequisite for an optimum filling result?

A root canal correctly prepared for obturation with warm gutta-percha should exhibit a continuously conical cavity – with the lowest diameter in the foramic area and the highest diameter in the coronal area. A higher conicity makes thermoplastic filling methods, such as the carrier-based technology, easier. A preparation with a conicity of 5% or higher meets these requirements.

11. Which conicity does the obturator have?

In the cold state the obturator has a conicity of 2–2,5%, depending on its size. Once the obturator is heated and inserted into the root canal, the gutta-percha is deformed and adjusts itself to the canal anatomy and its respective conicity.

FAQ — Frequently Asked Questions

12. Can GuttaMaster® also be used for warm obturation in curved canals?

Yes. GuttaMaster® is a very effective filling method for curved root canals.

Especially curved canals often pose the problem that with normal warm gutta-percha filling the used heat-conductors and pluggers often cannot be inserted to full working length due to the challenging anatomy. This results in the gutta-percha being not warm enough in the apical region and an insufficient plastification of the root filling.

The GuttaMaster® system however allows for a warm gutta-percha filling up to full working length, even in curved canals. The high flexibility of the GuttaMaster® carrier enables an application of gutta-percha also in strongly curved canals without exposing the plastic core.

13. What is the best way to fill two communicating root canals?

Option A: Simultaneously insert two GuttaMaster® obturators into the two root canals.

NOTE: Two obturators cannot be heated simultaneously. However, you can leave the first obturator in the oven until the second obturator is heated up. In the meantime the optimal heating temperature for the obturator in the first holder will be maintained for 90 seconds in the meantime.

Option B: Place a verifier in one canal as a placeholder, thus minimizing the transfer of gutta-percha when filling the other root canal. After having filled the first canal, remove the verifier and fill the second canal with another GuttaMaster® obturator.

14. Is the carrier coated with gutta-percha normally used in dental medicine?

Yes. GuttaMaster® obturators are coated with alpha gutta-percha. A special heat treatment lowers the melting point of the gutta-percha, enhancing its ease of flow. It is therefore better suitable for warm obturation methods than beta gutta-percha.

15. Can it happen that no gutta-percha flows into the apical area?

If you stick to the sequence of steps described in the user tips, both the apical area and the lateral canals can be filled with gutta-percha.

Kindly note that this filling method should be trained, before using it on patients, to enhance confidence when applying the GuttaMaster® technology. If you are interested in training options, please contact your local distributor.

16. Can it happen that the gutta-percha is stripped off the carrier and the carrier then has direct contact with the root canal wall?

Yes, that is possible, for example if the chosen obturator was not heated optimally. Or if gutta-percha was accidentally wiped off when the obturator was lifted out of the oven and this obturator was then used for filling.

Take the GuttaMaster® obturator carefully out of the holder of the GuttaMaster® oven. Take care that the gutta-percha is not accidentally wiped off the plastic carrier and insert the obturator into the root canal with a slow, steady, apically-directed movement. It is important that the direction of movement always runs along the longitudinal axis of the carrier to prevent it from bending. If you follow these steps, you can prevent the gutta-percha from being whipped off and the entire root canal can be filled tightly.

Furthermore we recommend using a surgical loupe or a dental microscope.

17. How much time do I have to bring the heated-up GuttaMaster® obturator into the root canal?

You have approx. 10 seconds (transfer time) before the obturator has cooled off to a temperature where the material can possibly no longer be inserted to the desired depth. Therefore there is no need to rush the insertion of the heated obturator.

If, however, you need longer or you have to take a short break, the GuttaMaster® oven is designed to maintain the respective heating temperature for 90 seconds.

An audible signal will remind you several times that there is still an obturator in the oven. As soon as this signal stops, you should fill the root canal with the obturator in the oven within 10 seconds. If you exceed the 10 second transfer time, you can let the obturator cool down and then heat it up again. However, the obturator must not be damaged and the gutta-percha coating must be intact.

FAQ — Frequently Asked Questions

18. How long does it take for the heated gutta-percha to cool off in the root canal?

The gutta-percha cools down within 2 to 4 minutes and becomes increasingly rigid.

During this time it is important that the handle is not moved. Instead, constant pressure should be applied for a few seconds towards the apex. If the obturator is pulled or turned, the filling result can be affected.

19. If the handle of the GuttaMaster® obturator is removed with a GuttaMaster® PREPI® Bur, can vibrations occur which might loosen the obturator in the root canal?

Before and while you separate the obturator handle with the appropriate drill, you should apply even pressure to the handle with a finger until it separates. This prevents the obturator from coming loose. As soon as the gutta-percha has cooled off sufficiently, a loosening of the filling is impossible.

20. What happens with any air which may exist in the root canal before and during the filling process?

Kindly take note that if the obturator is inserted too quickly and with too much pressure, the air can be compressed to such an extent that the patient feels pain. To prevent this, always insert the heated obturator steadily and slowly into the root canal.

In order to protect the patient from feeling any pain, we recommend administering an anesthetic before the endodontic procedure.

21. In a tooth with multiple roots, how can I prevent gutta-percha from accidentally flowing into another canal?

After having prepared the root canals, seal off those canals not to be filled in the first step, for example by using a foam pellet. This blocks a possible transfer of the sealer or gutta-percha into other prepared root canals.

22. Can I heat up two obturators in the GuttaMaster® oven at the same time?

If you want to heat up several obturators, you can use the right and left holders alternately for increased efficiency. Wait for the signal tone of the first holder before going on with the heating of the second holder.

NOTE: Two obturators cannot be heated simultaneously. However, you can leave the first obturator in the oven until the second obturator is heated up. In the meantime the optimal heating temperature for the obturator in the first holder will be maintained for 90 seconds.

The GuttaMaster® instruments

GuttaMaster® Starter Kit

The GuttaMaster® Starter Kit contains all you need to get started with the GuttaMaster® filling technology:

V04 0501 000 000	1	GuttaMaster® Starter Kit with 220 V GuttaMaster® Oven
V04 0501 000 530	1	GuttaMaster® Starter Kit with 110 V GuttaMaster® Oven
V04 0508 025 236	6 pcs.	GuttaMaster® Plastic Verifiers for Mtwo®, Selection: 1 x each size 20, 25, 30, 35, 40, 45
V04 0506 025 520	20 pcs.	GuttaMaster® Endodontic Obturators for Mtwo®, Selection: 3 x 20, 4 x 25, 4 x 30, 3 x 35, 3 x 40, 3 x 45
V01 0558 029 230	120 pcs.	Paper Points, Selection: Sizes 15-40, Taper .04
V01 0559 029 230	120 pcs.	Paper Points, Selection: Sizes 15-40, Taper .06
V04 1020	1	2Seal root canal sealer
V04 0515 000 000	6 pcs.	GuttaMaster® PREPI® Burs, Selection: 1 x \emptyset 0.70, 1 x 0.90, 2 x 1.2, 2 x 1.5
V04 0509 000 000	1	GuttaMaster® Practice Kit: 2 plastic practise blocks 2 obturators

GuttaMaster® Oven

V04 0500 000 000	1	GuttaMaster® Oven 220 V
V04 0500 000 530	1	GuttaMaster® Oven 110 V

GuttaMaster® Plastic Verifiers for Mtwo®

V04 0508 025 020-045		6-pack individual selection, Sizes 20, 25, 30, 35, 40, 45
V04 0508 025 236	6 pcs.	Selection: 1 x each size 20, 25, 30, 35, 40, 45

**GuttaMaster®
Endodontic Obturators for Mtwo®**

V04 0505 025 020-045		6-pack individual selection, Sizes 20, 25, 30, 35, 40, 45
V04 0505 025 236	6 pcs.	Selection: 1 x each size 20, 25, 30, 35, 40, 45
V04 0506 025 520	20 pcs.	Selection: 3 x size 20, 4 x 25, 4 x 30, 3 x 35, 3 x 40, 3 x 45

Paper points

V01 0558 029 230 120 pcs. Selection:
 Sizes 15-40, Taper .04

V01 0559 029 230 120 pcs. Selection:
 Sizes 15-40, Taper .06

Root canal sealer

V04 1020 1 2Seal, 4 ml each
 Paste A and B

GuttaMaster® PREPI® Burs

V04 0515 000 000 6 pcs. Selection:
 1 x \varnothing 0.70, 1 x 0.90,
 2 x 1.2, 2 x 1.5

Training material

V04 0509 000 000 1 GuttaMaster® Practice Kit:
 2 plastic practise blocks
 2 obturators

Special advantages of the GuttaMaster® filling technology:

- apical sealing of the canal including lateral canals and isthmi
- warm gutta-percha up to the end of root canal preparation
- simple filling of straight and curved root canals
- time-saving
- user-friendly

In the following cases an obturation with GuttaMaster® is not recommended:

- for teeth without sufficient apical constriction, e.g. with open apex, resorption, etc.
- with patients who might have an allergic reaction to sealer and/or gutta-percha
- cases in which a root canal preparation with a conicity of at least 5 % is not possible

VW000266 Rev. 0/29.05.09

VDW GmbH

P.O. Box 830954 • D-81709 München
Tel. +49 89 62734-0 • Fax +49 89 62734-304
www.vdw-dental.com • info@vdw-dental.com

Endo Easy Efficient®