


# Salli Saddle Chair dimensions

## Twin

Total weight 10,4 kg  
Seat weight 6,8 kg


Standard gas spring:

short: 465 ~ 600 mm  
medium: 590 ~ 770 mm  
long: 620 ~ 880 mm



Foot Height Control gas spring:

short: 455 ~ 590 mm  
medium: 575 ~ 765 mm  
long: 630 ~ 890 mm




User height (cm)	<160	160-180	>165
Gas spring	short	medium	long

## Swing

Total weight 10,2 kg  
Seat weight 6,6 kg


Standard gas spring:

short: 490 ~ 625 mm  
medium: 600 ~ 785 mm  
long: 645 ~ 905 mm



Foot Height Control gas spring:

short: 480 ~ 615 mm  
medium: 585 ~ 775 mm  
long: 655 ~ 915 mm




User height (cm)	<165	165-185	>170
Gas spring	short	medium	long

## MultiAdjuster

Total weight 11,1 kg  
Seat weight 7,5 kg


Standard gas spring:

short: 480 ~ 615 mm  
medium: 585 ~ 770 mm  
long: 635 ~ 895 mm



Foot Height Control gas spring:

short: 465 ~ 600 mm  
medium: 575 ~ 760 mm  
long: 640 ~ 900 mm




User height (cm)	<160	160-180	>165
Gas spring	short	medium	long

## SwingFit

Total weight 11,0 kg  
Seat weight 7,4 kg


Standard gas spring:

short: 500 ~ 635 mm  
medium: 610 ~ 800 mm  
long: 665 ~ 925 mm



Foot Height Control gas spring:

short: 490 ~ 625 mm  
medium: 605 ~ 790 mm  
long: 665 ~ 925 mm




User height (cm)	<165	165-185	>170
Gas spring	short	medium	long

## Professional

Total weight 9,7 kg  
Seat weight 6,0 kg


Standard gas spring:

short: 490 ~ 625 mm  
medium: 600 ~ 780 mm  
long: 645 ~ 905 mm



Foot Height Control gas spring:

not available  
not available  
not available




User height (cm)	<160	160-180	>165
Gas spring	short	medium	long

## Job

Total weight 6,0 kg


Standard gas spring:

short: 435 ~ 565 mm  
medium: 510 ~ 705 mm  
long: 590 ~ 845 mm



Foot Height Control gas spring:


not available  
not available  
not available



User height (cm)	<150	150-170	>160
Gas spring	short	medium	long


# Salli Saddle Chair dimensions

## Twin & Allround




User height (cm)	<165	165-185	>170
Gas spring	short	medium	long

## Multijuster & Allround




User height (cm)	<163	165-185	>170
Gas spring	short	medium	long

## Swing, Allround & big base




User height (cm)	<170	170-190	>175
Gas spring	short	medium	long

## SwingFit, Allround & big base




User height (cm)	<170	170-190	>175
Gas spring	short	medium	long

## Twin, Stretching support & big base



User height (cm)	<165	165-185	>170
Gas spring	short	medium	long

## Multijuster, Stretching support & big base




User height (cm)	<165	165-185	>170
Gas spring	short	medium	long

# Salli Saddle Chair dimensions



## Classic

Total weight 11,0 kg  
Seat weight 7,4 kg

### Standard gas spring:




### Foot Height Control gas spring:




User height (cm)	<160	160-180	>165
Gas spring	short	medium	long

## Classic & Allround

### Standard gas spring:




### Foot Height Control gas spring:




User height (cm)	<165	165-185	>170
Gas spring	short	medium	long

## Classic & Stretching Support & big base

### Standard gas spring:



### Foot Height Control gas spring:



User height (cm)	<165	165-185	>170
Gas spring	short	medium	long

## Accessories - further information

Accessories	weight	notify
Allround	4,36 kg	big base with swing and swingfit
Elbow rest	2,54 kg	
Stretching Support	2,16 kg	big base
Surgeon	4,32 kg	
Ergorest with stretching support	4,00 kg	
Foot adjustment	1,28 kg	
Gas spring lengths (cylinder)	length	notify
mini	185 mm	available on request only
short	235 mm	
medium	325 mm	
long	385 mm	